

Amphoteric Compound Chemical Eye/Skin Splash Decontamination: Clinical Experience

Hall AH¹, Mathieu L², Burgher F², <u>Fosse C²</u>, Blomet J²

¹Toxicology Consulting and Medical Translating Services, Inc., Laramie, Wyoming, USA and Colorado School of Public Health, Denver, Colorado, USA ²Laboratoire Prevor, Valmondois, France



Objective and method

- To evaluate clinical experience with using an amphoteric compound in a hypertonic solution as a decontamination solution for eye/skin chemical splashes, extending previously published data.(1-4)
- A review of cases of clinical use of an amphoteric solution for eye/skin chemical splash decontamination voluntarily reported to the manufacturer as part of a post-marketing surveillance program was done.

Hall AH et al. Vet Human Toxicol 2002;44:228-231.
Merle H et al. Burns 2005;31:205-211.
Nehles J et al. Cutan Ocular Toxicol 2006;25:249-258.
Donoghue AM. Int J Dermatol 2010;49:894-900.



Mechanism of the chemical injury

- = result of a contact of the skin/the eye with an irritant or a corrosive
- The severity of the burn depends mainly on:
 - nature and concentration of the chemical involved,
 - energy level involved,
 - time of contact.



Caustic soda (NaOH) burn at 24 hr not washed

- And also:
 - physical factors like the pressure or the temperature
 - the Total Body Surface Area (TBSA)
 - and whether the skin/eye is healthy or not.
- Chemicals can also be toxic! (Ex: hydrofluoric acid, HF)



What is Diphoterine[®]?

- An active washing solution for ocular and cutaneous chemical splashes
- An <u>aqueous solution</u> containing the fundamental properties of water:
 - mechanical effect of pulling the chemical agent away from the surface of the body tissue.

an amphoteric solution

 Acts in the same way on the acids and the bases, with a rapid return towards a zone of physiological pH. Also amphoteric for oxidisers/reducing agent and solvents.

A hypertonic solution

 Stops the penetration of corrosive chemicals into the tissues creating a flux from the inside to the outside of the body.

A medical device EC 0459, Ila class



Collected results

> From 1987 to 2010

- 58 reports from 52 different industrial facilities/organizations
- >44 total cases reported
- Some cases exposed several people on several body areas





Chemicals involved

No.	Chemical	Range of conc.	Ν	Conc. not specified
	Sodium hydroxide (NaOH)	flakes – 50% - 25%	14	4
-	sulfuric acid (H ₂ SO ₄)	98% - 96% - 87%	9	3
	Caustic solution or basic substance	pH>8.2	6	5
1	Phenol & derivates		4	4
	Nitric acid (HNO ₃)	100% - 70%	2	
	Lubricants		2	2
	Trichloroethane		1	1
-	Acrylic acid	100%	1	
	Acrylamide	50%	1	
	Acrylic varnish		1	1
(L	DMEA		1	1
-	Others		2	1
	Total		44	22



Body surfaces splashed

- 1 eye: 8 2 eyes:4 **16 eyes**
- **Face: 11**
- Forearms: 3
- Hands: 3
- Arms: 3
- Legs: 2
- Feet: 2
- Neck: 2
- Cheek: 1
- Chest: 1
- Back: 1
- Total body: 1
- Unknown: 4

- Eyes and face were the most involved body areas
- Skin splashes(5): Mean of TBSA: 9.1%

FIRST AID MANAGEMENT

- 40 cases washed with Diphoterine[®] first
- 3 cases washed with tap water first
 - 1 phenol splash washed with tap water + PEG + Diphoterine[®]

PEG: Polyethylene Glycol

⁵ Management of the major burn, Ames WA, practical procedures, Issue 10(1999) art. 10

	First alu management				
4		Yes	No	Not mentioned	
	Secondary treatment	4	9	31	
2	Lost work time	1	26	17	
1	Sequelae	9	19	16	



- In 2 cases, Diphoterine[®] was only used on *some but not all* exposed body parts. Outcome was worse in those parts not washed with Diphoterine[®].
- In one case, Diphoterine[®] was used rapidly on some body parts and only erythema developed. However, sodium hydroxide soaked through a shoe and Diphoterine[®] washing was delayed by 15 minutes, resulting in a second degree skin injury (see picture).
- In the 3 cases with water washing first, the outcome in 2 of 3 cases was generally less good than in the cases where Diphoterine[®] was used first.



Other data(4)

- Clinical case series of 180 alkali splashes on skin from Oct. 2006 to March 2008
- Dr M. Donoghue, Chief Medical Officer, Alcoa Australia

	Diphoterine [®] first	Water first
n	138	42
Time elapsed	1 min.	5 min.
No chemical burn	52.9%	21.4%
Blisters or more severe signs	7.9%	23.8%

4 Donoghue AM. Int J Dermatol 2010;49:894-900



Conclusion

Diphoterine[®] is currently used in industry, as a first aid solution

Used as primary action, it limits chemical burn development

In this non-random retrospective series, Diphoterine[®] was associated with good clinical outcomes.



Thank you

