

# A Review about Diphotérine®

## The solution for emergency decontamination of eye/skin chemical splashes

Presented at the Occupational Hygiene 2000, Manchester, UK, 11-13 April 2000

Mathieu, L., Blomet, J., PREVOR Laboratory, Valmondois, France (<http://www.prevor.com>); Girard, M., Rhône-Poulenc, La Rochelle, France ; Uellner, H., Martinswerk, Bergheim, Germany ; Nehles, J., Mannesmann, Remscheid, Germany

### Introduction

Diphotérine® is an eye/skin decontamination solution for chemical splashes<sup>(1)</sup>, produced by the PREVOR laboratory<sup>(2)</sup>, France. Its chemical and physical properties allow a quasi polyvalent rinsing of chemical splashes with a quick return to a physiological state. Diphotérine® stops the penetration of the chemical product<sup>(3),(4)</sup>. The emergent use of Diphotérine® will avoid any damage.

### Materials and methods

Most of the companies, mentioned in this report, were previously using water for emergency first aid decontamination of eye/skin chemical splashes without complete success : irremediable sequelae, numerous secondary care and loss of work. Subsequently, the Medical and Health and Safety Services decided to introduce Diphotérine® (or Previn®) for rinsing chemical splashes and to train workers to use it correctly. Previn® is the German version of Diphotérine®. Each ocular or cutaneous chemical splash was rinsed in emergency (some seconds to a few minutes) with Diphotérine®, on location, while undressing if necessary. Then each person went to the nursery for an examination. In the Mannesmann company, a secondary rinsing with Diphotérine® was performed in the nursery.

For eye splashes, Sterilized Individual Eye wash (SIEW) of 50 mL was used within 10 seconds while 500 mL of Diphotérine® was used for a longer time of contact (about 1 minute). For small skin splashes such as a hand or a face, a spray of 100 or 200 mL of Diphotérine® was used depending of the time of contact and the area. A big skin splash such as a leg or the chest involved the use of an autonomous portable shower with 5 L of Diphotérine® in the minute following the accident.


### Results

#### TOXICITY

Diphotérine® is a non toxic solution (Oral Toxicity LD<sub>50</sub>>2000 mg/Kg, Test, 6564TAR 1990 CIT; Acute Dermal Toxicity LD<sub>50</sub>>2000 mg/Kg, test 133/9, 1988, Safepharma). Its residues with acids and bases are non irritant (tests 6463TAL/6462TAL, 1990, CIT)<sup>(5)</sup>.

#### IN VITRO RESULTS

An in vitro experimental protocol<sup>(6),(7)</sup> has been adjusted by the PREVOR Laboratory for simulating ocular external rinsings of chemical splashes. It shows the greater efficacy of Diphotérine® rinsing versus water or saline solution rinsing of chemical splashes.



### CASE REPORTS

#### Corrosive splashes rinsed with diphotérine®

year	No cases	Firm /country	Exposure	Body Surface Area
1999	1	Knoll AG Germany	96% sulfuric acid	cheek
1998	1	Giesecke&Debrient Germany	100% nitric acid	hand
1995	1	Metaleurop Germany	96% sulfuric acid	Face + neck
1993	1	Mewa Germany	50% soda*	forearm
1991	1	Aluisuisse France	Soda flakes	Left eye
1991	2	Orgachim France	98% sulfuric acid	Face+neck+shoulders, (id + legs)

\* : soda = sodium hydroxide, \*\* : Cream ointment for the first exposure


**Result : no sequelae, no secondary care\*\*, no loss of work**

#### A SERIE<sup>(8)</sup> OF 24 CHEMICAL SPLASHES RINSED WITH DIPHOTÉRINE® in the Mannesmann factory Remscheid, Germany (1994-1998)

No Cases	Exposure	Body Surface
2	Quicklime (Calcium oxide)	One eye
1	30% soda (sodium hydroxide)	One eye
1	30% basic solution	One eye
6*	20% sulfuric acid	One eye
2*	5% phosphoric acid/35% nitric acid	One eye
2	Pure Sulfamic acid	One eye
1	5% sulfuric acid/ 7% nitric acid	One eye
4	20% sulfuric acid	Check/thorax/face/hand
1	53% nitric acid	head
1	75% phosphoric acid	Thorax+Genitals+right thigh
2	15% phosphoric acid	Forearm/hand
1	45% soda (sodium hydroxide)	Knee

**Results : no damage, no secondary care, no loss of work except 2 accidents\* with one day lost from work**

#### A STATISTICAL STUDY<sup>(7)</sup> ABOUT 180 CHEMICAL ACCIDENTS in the Rhône-Poulenc factory, La Rochelle, France (1987-1992)



**Results : using Diphotérine® instead of water completely suppressed loss of work and significantly reduced the need for secondary care**

#### A STATISTICAL STUDY<sup>(7)</sup> ABOUT 42 SODA (40G/L TO 600 G/L SODIUM HYDROXIDE) SPLASHES RINSED WITH DIFFERENT WASHING SOLUTIONS in the Martinswerk factory, Bergheim, Germany (1991-1993).

	Diphotérine®	Acetic acid	Water
No days lost from work	0,18d ± 0,4	2,91 d ± 4,3	8 d ± 8,12
No care	100% ± 15%	0% ± 15%	0 ± 15%
Simple care	0 ± 15%	80% ± 15%	25% ± 15%
Medical care	0 ± 15%	20% ± 15%	75% ± 15%

**Results : Using Diphotérine® resulted in a noticeable decrease in sick leave average and a standard deviation. No secondary care was necessary. There is a significant difference (p<0.05) between Diphotérine® and water concerning secondary care and loss of work.**

### Conclusion

The emergency use of Diphotérine® is a good way for the decontamination of ocular or cutaneous chemical splashes. It achieves a reduction of loss of work and secondary care in all cases and avoids the sequelae for the workers.

### References

- (1) Falcy M, Blomet J, DMT 1997,70, 137-146
- (2) Burgher F, Blomet J, Mathieu L 1996 Le Risque Chimique et la Santé au Travail, Ed PREVOR, France, ISBN2-9510211-0-0
- (3) Gérard M, Louis V, Merle H, Josset P, Menerath JM, Blomet J, J. Fr. Opht 1999, 22, 1047-1053
- (4) Josset P, Meyer MC, Blomet J, SMT, 1986, 25-33
- (5) Langefeld S, Blomet J, Mathieu L, Schrage N, Kompa S, Tympner J, Toxicology Letters 1999, 109/Suppl 1, 97-98
- (6) Josset P, Pelosse B, Saroux H, Bull Soc Opht Fr 1986, 6-7, 765-769
- (7) Oral presentation First International Congress Evolution of the Knowledge of Chemical Burns 1997, Oct 16-17, La Baule, France
- (8) Mathieu L, Blomet J, Nehles J, Dupont V SFUM/STC, 1999 April 22, Marseille, France