

Diphoterine® eye/skin chemical splash decontamination solution:

An updated review of safety and efficacy data

Hall AH¹, Mathieu L², Blomet J²

¹Toxicology Consulting and Medical Translating Services, Azle and Springtown, Texas and Colorado School of Public Health, University of Colorado-Denver, Denver, Colorado, USA; ²Laboratoire PREVOR, Valmondois, France

Accepted for poster presentation at the APAMT 2016 Congress, Singapore, 2016.

Innocuity of Diphoterine® solution


Test	Results
Ocular irritation	Non irritant
In vitro Evaluation of the eye	No cytotoxic or irritant potential
Cutaneous irritation	Non irritant
Ocular irritation of a residue	Non irritant
Ocular irritation of a residue	Non irritant
Oral toxicity	Non toxic, LD ₅₀ > 2000 mg/kg
Acute dermal Toxicity	Non toxic, LD ₅₀ > 2000 mg/kg
Sensitisation	Non sensitising
Mutagenesis	Non mutagenic
Cytotoxicity	Non cytotoxic
Anti-inflammatory potential	Non anti-inflammatory
Local tolerance on damaged skin/healthy skin	No irritant or toxic effects
Local skin tolerance	Non irritant

Formulation, classification of Diphoterine® solution

Aqueous saline solution containing Diphoterine®, amphoteric agent	
Medical device, Class IIa, CE 0459	can be used on healthy and damaged tissues
Does not contain phosphates	pH ranging between 7.2 and 7.7
Limpid and colourless liquid	Density : 1.034
Osmotic pressure : 820 mosmoles/kg	Sterile solution (by autoclave)

In vivo study in the rabbit eye:

Burn due to concentrated ammonia



No stromal edema and significant decrease of pH with an inflection of the curve when washing with Diphoterine® solution

An In Vivo cutaneous prospective, randomized, blind study:

Burn due to concentrated hydrochloric acid

Biomarkers	Results
IL ₆	Significantly decreased by the washing with Diphoterine® solution compared to other washing solutions at 48 hours and 7 days, 0.001 < p < 0.05
Substance P	Decreased and significantly different with Diphoterine® solution versus other groups at 6 and 48 hours, p < 0.05
B-Endorphin	Significantly increased when washing performed with Diphoterine® solution compared to no washing or other washing solutions after 7 days, p < 0.05

The differences in serum levels of biological markers and wound healing were likely due to the superior washing properties of Diphoterine® solution such that less HCl was left on the skin to produce injury (Cavallini, 2004).

Case Reports: Summary of observations when washing with Diphoterine® solution during an industrial chemical exposure (2014)

53 firms gave a testimony	Firms have described the following effects:
64 victims of specified accidents	- Decrease in secondary care
8 testimonies of « Everyday » use without special accident stories	- Improvement of the symptoms during washing
16 accidents involving Acids	- Less pain
25 accidents involving Alkalis	- Decrease in severity of the burn injuries
37 eye exposures	- Decrease in sequelae
35 skin exposures	- Decrease in lost Time of Work
8 combined eye and skin splashes	

A retrospective comparative study at the hospital (Verbelen 2016): Rinsing with Diphoterine® or Hexafluorine® versus water

Title / Author	Comparison	Dipho/Hexa	Water	Results
Chemical injury: 4 years of experience with an advanced approach	Need for surgery (number)	5	43	p < 0.0001
Verbelen	Average hospital stay per patient (Days)	3.48	7.76	P = 0.031
112 patients				

Less need for surgery, shorter hospital stay, patient go back earlier to work

Clinical data

CUTANEOUS BURN STUDIES											
Title Author	Study Type	Chemical	Time to rinse (minutes)	Area burned (%)	Burn outcome	Intervention Outcome	Observations / Conclusions				
The clinical efficacy of Diphoterine® in the management of cutaneous chemical burns: A 2-year evaluation study Zack-Williams 131 patients	Comparative Study	"Alkali burns"	Diphoterine (Dipho)	Water (H ₂ O)	Dipho H ₂ O	Dipho H ₂ O	The Diphoterine-treated patients were significantly younger (37.7 vs 43.2 years, p = 0.044) than those treated without. Patients who received Diphoterine presented to hospital significantly earlier than those who did not receive it (0.5 vs 2.55 days p = 0.006). There was a significant change in the wound pH pre- and post-application of Diphoterine® solution, compared to patients who were treated with water irrigation only, with a pH change of 1.076 vs 0.4 (p < 0.05).				
		n = 72 cases (55.0%)	0.57	2.15	1.76	1.25		8.07	7.77		
		"Acids"	n = 24 cases (18.3%)	n = 47	n = 84						
		"Other chemicals"	n = 35 cases (26.7%)	p = 0.004*		p = 0.203			p = 0.369		
Diphoterine® for alkali chemical splashes to the skin at alumina refineries Donoghue 180 patients	Comparative Study	"...strong alkali solutions (primarily sodium hydroxide)"	Dipho first	Water first, then Dipho	Dipho H ₂ O	Severity was recorded by medical personnel in the initial assessment. These results were not published with the article.	Severity	Dipho	H ₂ O		
		n = 135	n = 42	2.9	11	1.6	2.9	1 (no sign)	73 cases	9 cases	
Diphoterine® for Emergent Decontamination of Skin/Eye Chemical Splashes: 24 Cases Nehles	Case Series	ACIDS	"...nearly immediate (within the first 30-120 seconds after exposure)"								
		HNO ₃ (53%) H ₂ SO ₄ (20%) H ₂ SO ₄ (20%) H ₃ PO ₄ (16%) H ₂ PO ₄ (15%) H ₂ PO ₄ (75%) H ₂ SO ₄ (20%)	Head Right cheek Throat Left forearm Face Right hand Thorax, genitals Right hand	Lost work days	0 0 0 0 0 0 0	No sequelae in any case.	No sequelae, no need for secondary care, no loss of work				
		BASES	NaOH (45%)	(same as acids)	Knee						

* statistically significant
** approaching statistical significance
TBSA = Total Body Surface Area

OCULAR BURN STUDIES

Title Author	Study Type	Chemical	Time to rinse (minutes unless otherwise listed)	Area burned	Burn outcome	Intervention Outcome	Notes
Diphoterine for emergent decontamination of skin/eye chemical splashes: 24 cases. Nehles	Case Study	H ₃ PO ₄ /HNO ₃ (5/30-35%)			Left eye (L)	0	
		H ₂ SO ₄ (20%)			Right eye (R)	0	
		NH ₂ SO ₄ H (Powder)			R	0	
		H ₂ SO ₄ (20%)			Not reported	0	
		NH ₂ SO ₄ H (Powder)			Not reported	0	
		H ₂ SO ₄ (20%)			R	1	
		H ₃ PO ₄ /HNO ₃ (5/35%)			L	1	
		H ₂ SO ₄ (20%)			L	0	
		H ₂ SO ₄ /HNO ₃ (5/35%)			L	1	
		H ₂ SO ₄ (20%)			R	0	
		H ₂ SO ₄ (20%)			L	0	
		H ₂ SO ₄ (20%)			L	0	
		BASES					
		NaOH (30%)			R	0	
		"Basic Solution" (30%)			R	0	
		Quicklime (CaO)			R	0	
		Quicklime (CaO)			L	0	

Title / Author	Chemical	Time to rinse (minutes)	Eyes	Severity	Time to reepithelialization (days)	Observations / Conclusions	
Martinière (French West Indies) Evaluation of the use of an amphoteric solution as the rinsing product. Merle H Comparative study during 4 years 66 cases	Alkali n = 32 (48.5%) Dipho Phys. p < 0.0001* Javel ¹ n = 10 (15.1%) Dipho Phys. No p-value reported "Other" ² n = 24 (36.4%) Dipho Phys. 20 (55.6%) 4 (13.3%)	Roper-Hall Modification of the Hughes Classification Grade I Grade II Grade III	Dipho Phys.	Dipho Phys.	Grade I Corneal epithelial damage	DAP Phys. 1.9 ± 1 11.1 ± 1.4	Injuries were first irrigated with their respective solution at an average of the times listed (units in minutes). A second irrigation occurred 5 hours after the accident (5.1 ± 4.3 h). Significant difference for grade I and II. Not enough grade III to compare. All Grade IV injuries were rinsed with the physiologic solution and were therefore not included in this table.
			p = 0.49, NS	Specific eye not recorded	Grade II Corneal haze, iris details visible	p = 10 ⁻⁷ , significant	
			p = 0.79, NS	Specific eye not recorded	Grade III Total epithelial loss, stromal haze, iris details obscured	p = 0.02*, significant	
			p = 0.64, NS	Specific eye not recorded	Grade III 1/3-1/2 limbal ischaemia	p = 0.21, NS	

Title Author	Study Type	Chemical	Experimental groups	Pain level (inside CS cloud)	Time Interval between CS exposure and arrival... ...at the 'ready-for-action' checkpoint	Residual Pain...	Observations / Conclusions
Prevention of CS "Tear Gas" eye and skin effects and active decontamination with diphoterine Bvrar 22 participants	Comparative Study	CS group	Exposed to only CS.	9.7±0.5	2:28±0:25	2.3±0.5	Pain was scored according to a 10-point scale. A control group with 200ml low-pressure spray containers filled with water were prepared; however, the but officers refused to use water sprays due to their previous bad experiences with water decontamination after CS exposure. An aqueous, hypertonic, amphoteric, and chelating solution used prior to entering a riot reduces pain and recovery time after CS exposure. Moreover, in cases of CS exposures, decontamination with the aqueous, hypertonic, amphoteric, and chelating solution reduces facial pain.
		Pre-exposure group	Faces sprayed with Diphoterine (200ml) just before CS exposure.	5.6±1.1*	1:26±0:44*	1.1±0.4*	
		Post-exposure group	Faces sprayed with Diphoterine (200ml) immediately after CS exposure.	9.1±0.4	2:30±0:48	1.4±0.7*	

Title Author	Study Type	Chemical	Time to rinse (minutes unless otherwise listed)	Area burned	Burn outcome	Intervention Outcome	Observations / Conclusions
An amphoteric rinse used in the emergency treatment of a serious ocular burn Gerard 1 Grade IV case report	Case Study	Ammonia (Alcali®): 15.3%, pH: 12.8	"...1 h after the accident."		Grade IV Roper-Hall Classification • Cornea opaque, iris and pupil obscured • >1/2 limbal ischaemia • Poor prognosis Visual acuity 2/20	Time to Re-epithelialization (days) "Progression to healing [began at]" "Total re-epithelialization" 21 180	Rinsing was enhanced by instillation of local anaesthesia with oxybuprocaine eye drops. No need for surgery. Among the different rinsing solutions available, Diphoterine® seems to be valuable even after a longer delay of more than 10 min.

References:

- Hardwicke J, Hunter T, Staruch R, Molemen N. Burns. 2012;38(3):383-7.
- Hall AH, Blomet J, Mathieu L. Vet Hum Toxicol. 2002;44(4):228-31.
- Mathieu L, Burgher F, Blomet J. J Chem Health & Safety. July/August 2007, 32-39
- http://www.prevor.com/fr/liste-des-produits-testes-et-lavages-recommandes-pour-les-solutions-diphoterine-et-hexafluorine
- Hall AH, Cavallini M, Mathieu L, Malbach H. Cutan Ocul Toxicol. 2009;28(4):149-56.
- Donoghue AM. Int J Dermatol. 2010;49(8):894-900.
- Nehles J, Hall AH, Blomet J, Mathieu L. Cutan Ocul Toxicol. 2006;25(4):249-58.
- Rihawi S, Frenzt M, Schrage NF. Graefes Arch Clin Exp Ophthalmol. 2006;244(7):845-54.
- Bvrar M. Human and Experimental Toxicology 2015, 1-6
- Viala B, Blomet J, Mathieu L, Hall AH. J Emerg Med. 2005;29(1):5-8.
- Merle H, Donnio A, Ayebova L, Michel F, Thomas F, Ketterle J, et al. Burns. 2005;31(2):205-11.
- Gérard M, Merle H, Chiambaretto F, Rigal D, Schrage NF. Burns 2002; 28: 670-673
- Zack-Williams SD, Ahmad Z, Molemen NS. Ann Burns Fire Disasters. 2015 Mar 31; 28(1): 9-12
- Cavallini M, Puggioni V, Gazzola R. J of Plastic Dermatology. 2010;6(2):145-147
- Schrage NF, Rihawi R, Frenzt M, Reim M. Klin Monatsbl Augenheilkd 2004, 221: 1-98
- Cavallini M, Casati A. European Journal of Anaesthesiology 2004, 21, 389-392
- Cavallini M, de Brocard F, Corsi MM, Fassati LR, Baruffaldi Preis FW. Annals of burns and fire disasters 2004, XVII, 2, 1-5
- Malbach H, Hall AH. Chemical Skin Injury, Ed. Springer, 2014, ISBN 978-3-642-39778-3
- Schrage NF, Burgher F, Blomet J, Bodson L, Gerard M, Hall AH, Jasset P, Mathieu L, Merle H. Chemical Ocular Burns, Ed. Springer, 2011, ISBN 978-3-642-14549-0
- Atiyeh B, Barret JP, Dahai H, Dutelle F, Fowler Ann, Enoch S, Greenfield E, Magnette A, Rode H, Zhao-fan X. Wound International 2014
- Merle H, Gérard M, Schrage N. J Fr Ophthalmol. 2008 ;31(5) :1-12
- Deutsche Ophthalmologische Gesellschaft. Leitlinie Nr. 8. Verletzungen des Auges und seiner Anhangsgebilde. 2011 :1-8
- Lynn DD, Zukin LM, Delavalle R. Cutan Ocul Toxicol. 2016 Aug 31:1-7
- Verbelen J, Hoeksma H, Monstrey S. Nordic Burn Meeting, Uppsala, Sweden June 2016, oral presentation, to be published

