

CUTANEOUS LOCAL TOLERANCE OF DIPHOTERINE®

in human volunteers after an occlusive application during 48 hours

Mathieu L¹, Burgher F¹, Hall AH^{2,3}, Maibach HI⁴

Presented at the SOT congress, March 2008, Seattle, USA

Diphoterine® is an eye/skin decontamination solution for chemical splashes which contains an amphoteric compound. Its chemical and physical properties allow a quasi polyvalent rinsing of chemical splashes with a quick return to a physiological state. The local tolerance of Diphoterine® was evaluated on human volunteers after a single and occlusive application during 48 hours.

MATERIAL AND METHODS

The skin local tolerance was evaluated in 55 volunteers by IDEA (Dermatological Institute of Aquitaine), (study n°ID-07/0477 // 1.01-48H) Martillac, France). The investigation was carried out at Solar Test International (STI) in Romania, according to the internal standard protocol and ethical committee. A single application of 0.02ml of Diphoterine® (batch N°. D970112B) was applied on the outside surface of the arm and maintained for 48 hours in contact with the skin, with the help of an occlusive patch (Haye's Chambers). 57 volunteers were recruited, male and female, with normal skin, without any dermatological lesion on the experimental area, were included in the study. There were 8 volunteers aged 65 to 70 years old but this had no influence on the results. The clinical evaluation was made 30 minutes after the patch removal and takes in account the erythema, the papules, the vesicles and blisters. According to their intensity, the value of the irritation ranges from 0 to 4. The total sum of the scores, divided by the number of volunteers, defines the average irritation index which allows the product to be classified (Table 2).

non irritant	I.I.M. ≤ 0.20
slightly irritant	0.20 < I.I.M. ≤ 0.5
moderately irritant	0.50 < I.I.M. ≤ 2
very irritant	2 < I.I.M. ≤ 3
severely irritant	I.I.M. > 3

Table 2: Irritation index

RESULTS

Of the 57 recruited volunteers, 55 were included in the study. The average index of irritation I.I.M. for Diphoterine® is 0.00.

DISCUSSION

The results of this test show that Diphoterine® did not induce skin irritation. It is in accordance with the results obtained in a previous test where local tolerance was evaluated in the rabbit on both normal and scarified skin, by non-occlusive or semi-occlusive application. Semi-occluded application was used to mimic contact of wet clothes soaked with Diphoterine® and non-occluded application for Diphoterine® itself on skin as recommended protocol. No irritating or toxic effect was observed.

Skin sensitization study in the guinea pig has been also performed following the Magnusson and Kligman method and published. Other toxicological tests have also been performed in order to evaluate the toxicity of Diphoterine® by other routes of exposure. The acute toxicity of the product Diphoterine® was also evaluated in the rat by the oral route. The LD50 of the product Diphoterine® administered by the oral route in rats was found greater than 2000 mg/Kg, a dose at which no signs of toxicity were observed : no deaths, no change

of general behavior or bodyweight and no evidence of macroscopic abnormalities when necropsy was performed. The Ames test has also been performed and Diphoterine® did not show any mutagenic activity in the bacterial reverse mutation test with *Salmonella typhimurium* and *Escherichia coli*.

The conditions of exposure, used for the evaluation of local tolerance in human volunteers, are much more severe than what is required by the standard protocol of the use of Diphoterine® for decontamination of chemical splashes: a unique dose within one minute of contact, 500 ml for washing one eye and 5 liters during 5 minutes for the complete body.

All these toxicological data are in accordance with the post-marketing surveillance program of Diphoterine® use where no adverse effects of this eye/skin chemical splash decontamination solution have been noted.

These innocuousness data open the range for the use of different dermal applications of Diphoterine®, such as delayed management of chemical splashes in victims of chemical assaults for which more and more countries are now concerned by this problem.

A comparative clinical study between saline solution versus Diphoterine® on eye burns has already shown the interest of its delayed use: for grade 1 and grade 2 burns, reepithelialisation time is shorter when rinsing has been done with Diphoterine®: 1.9 ± 1 day, compared with 11.1 ± 1.4 days (p = 10-7) and 5.6 ± 4.9 days compared with 10 ± 9.2 days (p = 0.02).

An experimental and comparative study of a concentrated hydrochloric acid burn on rat skin showed that skin flushing with Diphoterine® stopped the burn evolution and reduced substance P release and higher concentrations of b-endorphin were observed compared to other washing solutions.

CONCLUSION Diphoterine® can be considered as non-irritant after a single application of 48 consecutive hours.

REFERENCES

- Mathieu L, Burgher F, Blomet J. *J Chem Health and Safety*, 2007, 14, 4, 32-39
- Mathieu L, Burgher F, Hall AH. *Cutan Ocu Toxicol* 2007; 26:181-187
- Merle H, Donnio A, Ayeboua L, Michel F, Thomas F, Ketterle J, Leonard C, Josset P, Gérard M. *Alkali ocular burns in Martinique (French West Indies): Evaluation of the use of an amphoteric solution as the rinsing product*. *Burns*, 2005, 31, 205-211
- Magnusson B, Kligman AM. *The identification of contact allergens by animal assay: the guinea pig maximization test*. *J Invest Dermatol* 1969, 52, 268-276
- Mathieu L, Burgher F, Hall AH. *Diphoterine® chemical splash decontamination solution: skin sensitization study in the guinea pig*. *Cutaneous and Ocular Toxicology*. 2007, 26, 181-187
- Milton R, Mathieu L, Hall A, Maibach HI. *Chemical Assault and Eye/Skin Burns: Two representative Cases*, Report from the Acid Survivors Foundation, and Literature Review, submitted for publication
- *Acute toxicity by oral route in the rat*, Study 6564 TAR, Centre International of Toxicology, Evreux, France, 1990
- *Bacterial reverse mutation test*, Study 29023 MMT, Centre International of Toxicology, Evreux, France, 2006
- Cavallini M, Casati A. *A prospective randomized, blind comparison of saline, calcium gluconate and Diphoterine® for washing skin acid injuries in rats: effects on substance P and b-endorphin release*. *European Journal of Anaesthesiology*, 2004, 21, 389-392

1 - PREVOR Laboratory, Valmondois, France.

2,3 - Toxicology Consulting and Medical Translating Services, Inc., Laramie Wyoming, USA; Department of Preventive Medicine, University of Colorado Health Sciences Center, Denver, Colorado USA.

4 - Howard I. Maibach, M.D., Department of Dermatology, University of California-San Francisco, San Francisco, CA, USA



PREVOR
PRÉVOIR ET SAUVER
Laboratoire de Toxicologie & Maîtrise du Risque Chimique