X



## **Technical data sheet**

# DIPHOTERINE® solution Wall-mounted Station - LMPE.DAS



## **Advantages**

- Fixed and easily identifiable.
- Easy and ready to use.
- Allows quick action without assistance.
- The entire face can be washed, with an LPMD eyewash for each eye and MINI DAP spray.
- Gentle washing without pressure. Autonomous system.
- Sterile solution.
- Easy control of the washing effectiveness: Use each single packaging in its entirety.
- Easy and adaptable installation.
- Watertight, in accordance with industrial hygiene and cleanliness.
- Without special care and without special installation (no need for electricity or evacuation, pipe or valve, no more leakage problems...).
- The washing station wall is therefore easily removable for independent work spot to modular nature.
- In conformity with EN 15454-4.
- EC certification.

### **Characteristics**

	/
Supplier	PREVOR Moulin de Verville 95760 Valmondois FRANCE +33-1-30-34-76-76
Sizes	Height: 420 mm Width: 330 mm.
Capacity	2 x 500 ml (LPMD) - DIPHOTERINE® solution eyewash 1 x 200 ml (LOA) - AFTERWASHII® solution eyewash 1 x 200 ml (MINI) - DIPHOTERINE® solution spray
Shelf-life	2 years
Content	DIPHOTERINE® solution, AFTERWASHII® solution.
Use	Up to 2 eyes and a face.
Maximum time of use recommended	1 minute maximum.

