

Evaluation of Lesions Caused by Hydrofluoric Acid (HF) on Human Skin *ex vivo* and Decontamination With Tap Water + Calcium Gluconate or Hexafluorine®

Mathieu L¹, Lati E², Burgher F¹, Gasser P², Hall AH³,
Peno-Mazzarino L², Maibach HI⁴, Blomet J¹

¹Laboratoire Prevor, Valmondois, France

²Center for Biological Research and Cutaneous Experiments (BIO-EC), Clamart, France

³Toxicology Consulting and Medical Translating Services, Inc., Laramie, Wyoming, USA and Colorado School of Public Health, Denver, Colorado, USA

⁴Department of Dermatology, University of California-San Francisco, San Francisco, California, USA


Summary

- Human skin explants *ex vivo* represent:
 - A means of testing hazardous substances and decontamination measures
 - *Without* using live animal models
 - In cases where studies in living humans are not
 - *Ethically feasible*
 - *Technically feasible*

Human Skin Explants *ex vivo*

- The BIO-EC Laboratory, Clamart, France has developed:
 - A technique for ethically using human skin explants
 - Maintained alive *ex vivo*
 - For the study of caustic or irritant substances and their decontamination

Human Skin Explants



Human Skin Explants *ex vivo*

Methods

■ Products

- HF: furnished by Study Sponsor; FLUKA, reference 47610; Lot No. 7125A; exact concentration = 73%
- Tap Water
- Calcium gluconate 2.5% gel; KAYS Medical, Lot No. 2001/09
- Hexafluorine®; Prevor Laboratory; Mini-DAP 200 mL; Lot No. F971201C; Expiration Date 04/03/2009

Human Skin Explants *ex vivo*

Methods – Control Explants

- 21 control Human skin explants, diameter approximately 10 mm, prepared from an abdominoplasty of a 45-year-old woman
- *After* obtaining informed consent
- Explants preserved alive in BEM (BIO-EC's Explant Medium)
 - At 37° C
 - Humidified atmosphere containing 5% CO₂

Human Skin Explants *ex vivo*

Methods – Control explants

- The 21 control explants were divided into 7 groups, 3 explants per group, with 70% HF exposure for:
 - 20 seconds
 - 1 minute
 - 2 minutes
 - 3 minutes
 - 4 minutes
 - 5 minutes

Human Skin Explants *ex vivo*

Methods

- 70% HF applied
 - Soaked into filter paper disk
 - Chosen to take up 30 μL HF

HF Application



Human Skin Explants *ex vivo*

Methods – Control Explants

- As compared with explants with *no* HF application:
 - **20 Seconds contact:** *no alterations* of epidermis or dermis
 - **1 minute contact:** *slight alterations* in the upper epidermal layer
 - **2 minutes contact:** *Clearly increased alterations* in the upper epidermal layer
 - **3 minutes contact:** *obvious lesions* in the epidermal layer
 - **4 minutes contact:** *very clear lesions* in the epidermal layer
 - **5 minutes contact:** *obvious lesions* in the epidermis and papillary dermis
- **In these experimental conditions, the initial cellular alterations are observed after 2 minutes of contact with HF**

Human Skin Explants *ex vivo*

Methods - Decontamination

- 86 Human Skin explants, approximately 10 mm diameter, obtained from an abdominoplasty of a 35-year-old woman
- *After* obtaining informed consent
- Preserved alive in BEM medium
 - At 37° C
 - Humidified atmosphere with 5% CO₂

Human Skin Explants *ex vivo*

Methods - Decontamination

- Explants divided into 5 groups:
 - **Group 1:** 20 explants; (*controls, no HF application*)
 - **Group 2:** 18 explants: 70% HF, 20 seconds contact (*untreated controls*)
 - **Group 3:** 16 explants 70% HF, 20 seconds contact (*Tap Water washing 15 minutes + Calcium Gluconate Gel 1 g/cm²*)
 - **Group 4:** 16 explants 70% HF, 20 seconds contact (*washing with 200 mL Hexafluorine® over approximately 5 minutes*)
 - **Group 5:** 16 explants 70% HF, 20 seconds contact (*washing with 2 x 200 mL Hexafluorine® -- 400 mL total -- over approximately 10 minutes*)

Human Skin Explants *ex vivo*


Methods - Decontamination

- Histopathology Sampling Times:
 - **Group 1:** T0, 20 seconds, 5 minutes, 10 minutes, 15 minutes, 30 minutes, 1 hour, 2 hours, 4 hours, 24 hours
 - **Group 2:** 20 seconds, 5 minutes, 10 minutes, 15 minutes, 30 minutes, 1 hour, 2 hours, 4 hours, 24 hours
 - **Group 3:** After end of washing: 5 minutes, 10 minutes, 15 minutes, 30 minutes, 1 hour, 2 hours, 4 hours, 24 hours
 - **Group 4:** After end of washing: 5 minutes, 10 minutes, 15 minutes, 30 minutes, 1 hour, 2 hours, 4 hours, 24 hours
 - **Group 5:** After end of washing: 5 minutes, 10 minutes, 15 minutes, 30 minutes, 1 hour, 2 hours, 4 hours, 24 hours

Tap Water Washing



Hexafluorine® Application



Human Skin Explants *ex vivo*

Results - Histopathology

■ ***Untreated explants:***

- Independent of the length of explant survival, no modifications of the epidermal or dermal structures observed.

Human Skin Explants *ex vivo*

Results - Histopathology

- ***70% HF-exposed explants, no decontamination:***
 - ***Clear alterations in the epidermis and papillary and inferior reticular dermis at 5 & 10 minutes***
 - ***Clearly increased alterations after 15 & 30 minutes and 1 hour***
 - ***Epidermal structures clearly necrotic at 24 hours***

Human Skin Explants *ex vivo*

Results - Histopathology

- **70% HF-exposed explants washed with Tap Water + Topical Calcium Gluconate Gel:**
 - After 5 & 10 minutes, ***no alterations*** of the epidermis or dermis
 - After 15 minutes, ***very clear cellular alterations*** in the epidermis and papillary and inferior reticular dermis
 - The alterations ***decrease*** after 30 minutes, but are ***more definitively observed*** after 1 & 2 hours
 - ***Slight edematous lesions*** observed in the epidermal basal layer after 4 hours which ***clearly increase*** after 24 hours

Human Skin Explants *ex vivo*

Results - Histopathology


- **70% HF-exposed explants washed with Hexafluorine® 200 mL:**
 - After 5 minutes, some clear epidermal and dermal alterations in a single explant; ***none*** in other explants
 - After 10, 15, & 30 minutes, ***no alterations*** in the epidermis or dermis
 - After 1 & 4 hours, very clear epidermal and dermal alterations in a single explant; ***none*** in other explants
 - After 2 & 24 hours, ***no alterations*** in the epidermis or dermis

Human Skin Explants *ex vivo*


Results - Histopathology

- ***70% HF-exposed explants washed with Hexafluorine® 200 mL x 2 (400 mL total):***
 - After 5, 10 & 30 minutes and 1, 2, 4, & 24 hours:
 - ***No alterations of the epidermal and epidermal structures were observed***


EXPLANT NON EXPOSE MORPHOLOGIE GENERALE
Non exposed explant general aspect



EPIDERME HF 5 MIN DE CONTACT
EPIDERMIS 5 MIN HF CONTACT




HF 24 H DE CONTACT NECROSE COMPLETE
HF 24 H CONTACT TOTAL NECROSIS



		Control group No exposure non washing (20 explants)	HF without washing (18 explants)	HF Water washing + Calcium Gluconate (16 explants)	HF Hexafluorine® washing 400 ml (16 explants)
T0	Epidermis	GM = Good Morphology			
	Papillary dermis	GM			
	Reticular dermis	BM			
20 s					
	Epidermis	GM			
	Papillary dermis	GM			
	Reticular dermis	GM			
Total time after exposure		20 s			
5 min	Epidermis	GM	PN + AC	GM	GM
	Papillary dermis	GM	PN + AC	GM	GM
	Reticular dermis	GM	PN + AC	GM	GM
Total time after exposure		5 min	5 min	15 + 5 = 20 min	10 + 5 = 15 min
10 min	Epidermis	GM	PN	GM	GM
	Papillary dermis	GM	AC	GM	GM
	Reticular dermis	GM	PN + AC	GM	GM
Total time after exposure		10 min	10 min	15 + 10 = 25 min	10 + 10 = 20 min
15 min	Epidermis	GM	PN + AC	PN + AC moderate	GM
	Papillary dermis	GM	PN + AC	PN + AC	GM
	Reticular dermis	GM	PN + AC	PN + AC	GM
Total time after exposure		15 min	15 min	15 + 15 = 30 min	10 + 15 = 25 min
30 min	Epidermis	GM	PN + AC	Some necrosed cells	GM
	Papillary dermis	GM	PN + AC	GM	GM
	Reticular dermis	GM	PN + AC	GM	GM
Total time after exposure		30 min	30 min	15 + 30 = 45 min	10 + 30 = 40 min
1 h	Epidermis	GM	PN + AC	GM	GM
	Papillary dermis	GM	PN + AC	GM	GM
	Reticular dermis	GM	PN + AC	GM	GM
Total time after exposure		1 h	1 h	1 h + 15 = 1 h 15	1 h + 10 = 1 h 10
2 h	Epidermis	GM	PN + AC	GM	GM
	Papillary dermis	GM	PN + AC	GM	GM
	Reticular dermis	GM	PN + AC	GM	GM
Total time after exposure		2 h	2 h	2 h + 15 = 2 h 15	2 h + 10 = 2 h 10
4 h	Epidermis	GM	PN + AC	Slightly edematous cells with mild acantholysis	GM
	Papillary dermis	GM	PN + AC	GM	GM
	Reticular dermis	GM	PN + AC	GM	GM
Total time after exposure		4 h	4 h	4 h + 15 = 4 h 15	4 h + 10 = 4 h 10
24 h	Epidermis	GM	Total necrosis	Very edematous cells with a very clear cytoplasm	GM
	Papillary dermis	GM	PN + AC	PN + AC	GM
	Reticular dermis	GM	PN + AC	Lesser alterations	GM
Total time after exposure		24 h	24 h	24 h + 15 = 24 h 15	24 h + 10 = 24 h 10


Morphology at 24h

Without washing




Total necrosis with grey cytoplasm, pyknotic nuclei and acidophilic cytoplasm in papillary and reticular dermis

Water + Ca Glu



Many edematous cells in basal epidermis with very clear cytoplasm and basal membrane disruption. Pyknotic nuclei and acidophilic cytoplasm in papillary dermis; same lesion but weaker in reticular dermis

Hexafluorine®



Normal morphology
in all layers

Human Skin Explants *ex vivo*

Conclusions

- In the experimental condition utilized, 70% HF penetrated in < 5 minutes and caused massive damage to various skin layers
- After 20 seconds of 70% HF contact, Tap Water washing + Topical Calcium Gluconate Gel allowed temporary slowing of HF action on tissues, but a new evolution of injury occurred after 24 hours

Human Skin Explants *ex vivo*

Conclusions

- Washing with 200 mL of Hexafluorine® over 5 minutes associated with better results at 5 minutes and 1 & 4 hours
- With an ***absence*** of HF injury at 4 & 24 hours
- Washing with 200 mL x 2 of Hexafluorine® (400 mL total) over 5 minutes was associated with ***no tissue injury*** after 5 minutes to 24 hours

Human Skin Explants *ex vivo*

Conclusions

- Hexafluorine® should be considered as the initial decontamination measure for HF skin splashes

