



Technical data sheet

DIPHOTERINE® solution Autonomous Portable Shower - DAPD



+33 1 30 34 76 76

Moulin de verville, 95760
Valmondois, France



Advantages

- Easily identifiable, the color of the DAP is green
- to make the difference with the fire extinguisher.
- To be used like an extinguisher
- as a consequence does not entail a specific training.
- Ready to use, permit a quick action
- A soft and micronised spray that wets a lot to permit a rinsing
- without trauma and a perfect distribution of the solution.
- The injured person can use the DAP alone.
- Long flexible that permits a complete rinsing of the body.
- Transportable. It permits the first aider to be protected as he will not be in contact with the splashed parts of the injured person.
- The system can be used when the DAP is laying.
- Easy control of the rinsing efficiency : complete and sole use.
- Sterile solution.
- The leading insures the integrity of the system.
- Does not require a specific use and maintenance (does not require a feed, piping or floodgate system, no more problem of leak).
- The DAP can always be transferred for the modulable work shops.
- Can be stocked in a protective box or in an anti freeze protective box.
- Refillable with the possibility of a maintenance contract.
- In conformity with EN 15154-3
- EC certification

Characteristics

Supplier	PREVOR Moulin de Verville 95760 Valmondois FRANCE + 33-1-30-34-76-76
Sizes	640 * 220 * 200 mm
Capacity	5 litres
Weight	11Kg
Propellant	Carbon dioxide (CO2)
Shelf-life	2 years
Content	DIPHOTERINE® solution
Use	The entire body surface
Delay to interfere	Maximum 1 minute

