



+33 1 30 34 76 76

Moulin de verville, 95760  
Valmondois, France



### Advantages

- Light and small dimensions.
- Easy to identify.
- Easy and ready to use.
- Permits a quick action.
- The ergonomic eyecup facilitates the eye opening.
- Soft and autonomous rinsing without pressure.
- Sterile solution.
- Sealed to assure the integrity of the system.
- Easy control of the rinsing efficiency : complete and sole use.
- Waterproof and compact.
- Does not require a specific installation and maintenance.
- In conformity with EN 15154-4.
- EC certification.

### Characteristics

<b>Supplier</b>	PREVOR Moulin de Verville 95760 Valmondois FRANCE +33-1-30-34-76-76
<b>Sizes</b>	Height : 124 mm Diameter : 64 mm
<b>Capacity</b>	50 ml
<b>Weight</b>	95g
<b>Shelf-life</b>	2 years
<b>Content</b>	DIPHOTERINE® solution
<b>Use</b>	One eye
<b>Delay to start the rinsing</b>	Maximum 10 seconds
<b>Available accessory</b>	Belt sheath

