



## DIPHOTERINE® solution MICRO-DAP spray



| Advantages                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                                                                                 |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>• Light and small dimensions.</li> <li>• Easy and ready to use.</li> <li>• Permits a quick action.</li> <li>• A soft and micronised spray that wets a lot to permit</li> <li>• a rinsing without trauma and a perfect distribution of the solution.</li> <li>• Sterile solution.</li> <li>• Easy control of the rinsing efficiency : complete and sole use.</li> <li>• Does not require a specific installation and maintenance.</li> <li>• In conformity with EN 15154-3</li> <li>• EC certification.</li> </ul> |                                                                                 |
| Characteristics                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |                                                                                 |
| <b>Supplier</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | PREVOR<br>Moulin de Verville<br>95760 Valmondois<br>FRANCE<br>+33-1-30-34-76-76 |
| <b>Sizes</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | Height :178 mm<br>Diameter : 40 mm                                              |
| <b>Capacity</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | 100 ml                                                                          |
| <b>Weight</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | 200g                                                                            |
| <b>Propellant</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | Nitrogen                                                                        |
| <b>Shelf-life</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 2 years                                                                         |
| <b>Content</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | DIPHOTERINE® solution                                                           |
| <b>Use</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | Hand, cheek ...                                                                 |
| <b>Maximum delay before rinsing</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | 1 minute maximum                                                                |



+33 1 30 34 76 76

Moulin de verville, 95760  
Valmondois, France

PREVOR

