



DIPHOTERINE® solution MINI-DAP spray



Advantages	
<ul style="list-style-type: none"> • Light and small dimensions. • Easy and ready to use. • Permits a quick action. • A soft and micronised spray that wets a lot to permit • a rinsing without trauma and a perfect distribution of the solution. • Sterile solution. • Easy control of the rinsing efficiency : complete and sole use. • Does not require a specific installation and maintenance. • In conformity with EN 15154-3 • EC certification. 	
Characteristics	
Supplier	PREVOR Moulin de Verville 95760 Valmondois FRANCE +33-1-30-34-76-76
Sizes	Height : 235 mm Diameter : 50 mm
Capacity	200 ml
Weight	290g
Propellant	Nitrogen
Shelf-life	2 years
Content	DIPHOTERINE® solution
Use	Hand and forearm, face, thigh ...
Delay to interfere	1 minute maximum

+33 1 30 34 76 76

Moulin de verville, 95760
Valmondois, France

PREVOR

