

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| PRODUCTS | CAS N° | TOXICITY | CHARACTERISTICS | DIPHOTERINE | HEXAFLUORINE | ANALYSIS REFERENCE |
|--|------------|--|---------------------|-------------|--------------|--------------------|
| (HEXA)FLUOROSILICIC ACID | 16961-83-4 | SEE HYDROFLUOSILICIC ACID | | | * | |
| 1,1',1"-NITRILOTRI(2-PROPANOL) | 122-20-3 | SEE TRIISOPROPANOLAMINE | | | | |
| 1,1,1,3,3,3-HEXAMETHYLDISILAZANE | 999-97-3 | CORROSIVE | | * | | 03,03,31,003/FDS |
| 1,1,13,3,3-HEXAFLUORO-2-PROPANOL | 920-66-1 | HARMFUL | | * | | 00,09,26,002 |
| 1,1,1-TRICHLOROETHANE | 71-55-6 | HARMFUL/DANGEROUS TO THE ENVIRONMENT | SOLVENT | * | | |
| 1,1,2,2-TETRABROMOMETHANE | 79-27-6 | TOXIC/IRRITANT | | * | | 03,11,04,002/MSDS |
| 1,1,2-TRICHLOROETHANE | 79-00-5 | HARMFUL | SOLVENT | * | | |
| 1,1,2-TRICHLOROTRIFLUOROETHANE | 76-13-1 | DANGEROUS TO THE ENVIRONMENT | | * | | |
| 1,1-DIMETHYLHYDRAZINE ANHYDROUS | 57-14-7 | FLAMMABLE TOXIC | | * | | 90,09,20,001 |
| 1,2,4-TRICHLOROENZENE | 120-82-1 | HARMFUL/IRRITANT | | * | | |
| 1,2-BROMOFLUOROBENZENE | 1072-85-1 | HARMFUL/IRRITANT/FLAMMABLE | | * | | |
| 1,2-DIBROMOPROPANE | 78-75-1 | HARMFUL | | * | | 95,05,30,010 |
| 1,2-DICHLOROETHANE | 107-06-2 | TOXIC/FLAMMABLE | SOLVENT | * | | 98,08,17,027 |
| 1,2-DICHLOROETHYLENE | 156-60-5 | FLAMMABLE, HARMFUL, Dangerous to the environment | SOLVENT | * | | 070104004-MSDS |
| 1,2-DICHLOROPROPANE | 78-87-5 | HARMFUL | SOLVENT | * | | 090204005/MSDS |
| 1,2-DIMETHOXYETHANE | 110-71-4 | HARMFUL | GLYCOL ETHER | * | | 98,08,17,039 |
| 1,2-EPOXYETHANE | 75-21-8 | SEE ETHYLENE OXIDE | | | | |
| 1,2-PROPANEDIOL | 57-55-6 | | ALCOHOL | * | | |
| 1,2-PROPYLENEDIAMINE | 78-90-0 | CORROSIVE/HARMFUL | BASE | * | | 03,03,18,001 |
| 1,3,5-TRICHLOROBENZENE | 108-70-3 | HARMFUL/IRRITANT | | * | | |
| 1,3,7-TRIMETHYLXANTHINE | 58-08-2 | SEE CAFFEINE | | | | |
| 1,3-BROMOFLUOROBENZENE | 1073-06-9 | HARMFUL/IRRITANT/FLAMMABLE | | * | | |
| 1,3-BUTANEDIOL | 107-88-0 | IRRITANT | SOLVENT | * | | |
| 1,3-CYCLOHEXANEDIONE | 504-02-9 | IRRITANT/HARMFUL | CETONIC SOLVENT | * | | 98,08,17,001 |
| 1,3-DI(TER-BUTYLPEROXYISOPROPYL)BENZENE | | IRRITANT | CROSS-LINKING AGENT | * | | 061226001/EXP |
| 1,3-DIAMINOPROPANE | 109-64-8 | IRRITANT/HARMFUL | SOLVENT | * | | 98,08,17,025 |
| 1,3-DIAMINOPROPANE | 109-76-2 | TOXIC/CORROSIVE | BASE | * | | 98,08,17,024 |
| 1,3-DICHLORO-2-PROPANOL | 96-23-1 | TOXIC | SOLVENT | * | | 98,08,17,026 |
| 1,3-DICHLOROBENZENE | 541-73-1 | HARMFUL, IRRITANT | SOLVENT | * | | 100412001/MSDS |
| 1,3-DICHLOROPROPANE | 142-28-9 | IRRITANT/HARMFUL/FLAMMABLE | SOLVENT | * | | 98,08,17,002 |
| 1,3-PROPANEDIOL | 504-63-2 | | ALCOHOL | * | | |
| 1,3-PROPANOLAMINE | 156-87-6 | CORROSIVE | | * | | 95,05,30,031 |
| 1,4-BROMOFLUOROBENZENE | 460-004 | HARMFUL/IRRITANT/FLAMMABLE | | * | | |
| 1,4-BUTANEDIOL | 110-63-4 | HARMFUL | SOLVENT | * | | |
| 1,4-DIAMINOBTUTANE | 110-60-1 | FLAMMABLE,CORROSIVE | BASE | * | | 91116001/MSDS |
| 1,6-DICHLOROHXANE | 2163-00-0 | | | * | | 95,05,30,002 |
| 1-AMINO-2-PROPANOL | 78-96-6 | CORROSIVE | | * | | 95,05,29,016 |
| 1-AMINONAPHTHALENE | 134-32-7 | SEE NAPHTHLAMINE | | | | |
| 1-CHLORO-2,3-EPOXYPROPANE | 106-89-8 | SEE EPICHLORHYDRINE | | | | |
| 1-CHLORO-2,4-DINITROBENZENE | 97-00-7 | TOXIC/DANGEROUS TO THE ENVIRONMENT | | * | | |
| 1-CHLORO-2,4-DINITROBENZENE | 97-00-7 | TOXIC/DANGEROUS TO THE ENVIRONMENT | | * | | 95,06,19,013 |
| 1-DECANOL | 112-30-1 | HARMFUL/IRRITANT | ALCOHOL | * | | |
| 1-DIMETHYLAMINO-2-PROPYLAMINE | 62689-51-4 | HARMFUL / CORROSIVE | | * | | 95,05,29,013 |
| 1-DODECANOL | 112-53-8 | IRRITANT | | * | | |
| 1-DODECENE | 112-41-4 | IRRITANT | | * | | 95,05,29,019 |
| 1-HEPTANOL | 111-70-6 | HARMFUL | SOLVENT | * | | |
| 1-HEXADECANOL | 36653-82-4 | IRRITANT | ALCOHOL | * | | |
| 1-HEXANOL | 111-27-3 | HARMFUL | SOLVENT | * | | |
| 1-METHOXY-2-PROPANOL | 107-98-2 | FLAMMABLE | | * | | 89,09,04,004 |
| 1-METHOXYPROPYLACETATE | 108-65-6 | IRRITANT | SOLVENT | * | | 01,07,12,001 |
| 1-OCTADECANOL | 112-92-5 | IRRITANT | | * | | 95,05,29,008 |
| 1-OCTANOL | 111-87-5 | IRRITANT | SOLVENT | * | | 070919005-MSDS |
| 1-PENTANOL | 71-41-0 | HARMFUL/FLAMMABLE | SOLVENT | * | | |
| 1-PROPANAMINIUM,N,N,N-TRIMETHYL-3-[(1-OXO-2-PROPENYL)AMINO]CHLORIDE, POL | 75150-29-7 | HARMFUL TO THE ENVIRONMENT | COPOLYMER | * | | 071130002-MSDS |
| 1-PROPANOL | 71-23-8 | FLAMMABLE | ALCOHOL | * | | |
| 2,2,2-TRIFLUOROETHANOL | 75-89-8 | HARMFUL/SUSPECTED CARCINOGEN/FLAMMABLE | | * | | 01,04,19,001 |
| 2,2'-AMINOETHOXYETHANOL | 929-06-6 | CORROSIVE / HARMFUL | | * | | 95,06,01,003 |
| 2,2-DIMETHYL-1,3-PROPANEDIOL | 126-30-7 | IRRITANT | DIALCOHOL | * | | 98,08,17,042 |
| 2,4,6-TRICHLOROPHENOL | 88-06-2 | HARMFUL/IRRITANT | | * | | 00,09,26,001 |
| 2,4,6-tris(DIMETHYLAMINOMETHYL)PHENOL | 90-72-2 | HARMFUL/IRRITANT | | * | | |
| 2,4-DI-TERT-BUTYLPHENOL | 96-76-4 | | | * | | 02,03,27,001 |
| 2,4-PENTADIONE | 123-54-6 | SEE ACETYLACETONE | | | | |
| 2,5-FURANDIONE | 108-31-6 | SEE MALEIC ANHYDRIDE | | | | |
| 2,6-DI-TERT-BUTYL-PARA-CRESOL | 128-37-0 | HARMFUL / IRRITANT | HARMFUL / IRRITANT | * | | 95,05,29,006 |
| 2,6-DI-TERT-BUTYL-PARA-CRESOL | 128-37-0 | HARMFUL/IRRITANT | | * | | 95,05,29,006 |
| 2,6-DI-TERT-BUTYLPHENOL | 128-39-2 | | | * | | 02,03,27,001 |
| 25% LITHIUM TRI-SEC-BUTYLHYDROBORATE IN THF | 38721-52-7 | CORROSIVE, Reacts with water | REDUCING AGENT | * | | 100616001/MSDS |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | | |
|--|------------|--------------------------------------|--------------------|---|--|-------------------|
| 2-AMINO-2-METHYL-1- PROPANOL | 124-68-5 | IRRITANT | | * | | 95,06,19,015 |
| 2-AMINO-2-METHYL-1,3- PROPANEDIOL | 115-69-5 | IRRITANT | | * | | |
| 2-BUTANONE | 78-93-3 | IRRITANT FLAMMABLE | | * | | 04,11,02,001/ECH |
| 2-BUTOXY-2-ETHOXYETHANOL | 112-34-5 | SEE DIETHYLENE GLYCOL MONOBUTYLETHER | | | | |
| 2-BUTOXYETHANOL | 111-76-2 | HARMFUL | SOLVENT | * | | |
| 2-CHLOROETHANOL | 107-07-3 | VERY TOXIC | | * | | 95,06,01,007 |
| 2-CHLOROETHYLECHLOROFORMIATE | 627-11-2 | TOXIC / IRRITANT | | * | | 95,06,02,004 |
| 2-CHLOROPROPANE | 75-29-6 | HARMFUL/FLAMMABLE | SOLVENT | * | | 070621001/MSDS |
| 2-CHLOROPROPANE | 75-29-6 | HARMFUL/FLAMMABLE | SOLVENT | * | | 070621001/MSDS |
| 2-DODECANOL | 10203-28-8 | IRRITANT | | * | | |
| 2-EHTG, dest, IBC WEISS KUPAL,950Kg | 7659-86-1 | SENSITIZER | | * | | 05,04,18,003/MSDS |
| 2-ETHOXYETHANOL | 110-80-5 | TOXIC | GLYCOL ETHER | * | | 98,08,17,034 |
| 2-ETHOXYPHENOL | 94-71-3 | HARMFUL/IRRITANT | | * | | 98,09,08,001 |
| 2-ETHYLHEXYL NITRATE | 27247-96-7 | HARMFUL | NITRATE | * | | 090204001/MSDS |
| 2-HEPTANOL | 543-49-7 | HARMFUL/FLAMMABLE | SOLVENT | * | | |
| 2-HYDROXY-2-METHYLPROPIONITRILE | 75-86-5 | TOXIC | DECOMPOSED | | | 100414001/MSDS |
| 2-METHOXY-1-METHYLETHYL ACETATE | 108-65-6 | IRRITANT | SOLVENT | * | | 03,03,31,002/FDS |
| 2-METHOXYETHYL ETHER | 111-96-6 | TOXIC | SOLVENT | * | | 98,08,17,030 |
| 2-METHOXYETHYLAMINE | 109-85-3 | CORROSIVE/FLAMMABLE | BASE | * | | 98,08,17,016 |
| 2-METHYL-1-PROPANOL | 78-83-1 | HARMFUL/FLAMMABLE | SOLVENT | * | | |
| 2-METHYL-2 PROPANOL | 75-65-0 | HARMFUL/FLAMMABLE | SOLVENT | * | | |
| 2-METHYL-2-BUTANOL | 75-85-4 | HARMFUL/FLAMMABLE | SOLVENT | * | | |
| 2-METHYLPROPANAL | 78-84-2 | HARMFUL/FLAMMABLE | REDUCING AGENT | * | | 070919004/MSDS |
| 2-NITRO-P-ANISIDINE | | SEE 4-METHOXY-2-NITROANILINE | | | | |
| 2-OCTANOL | 4128-31-8 | IRRITANT | ALCOHOL | * | | |
| 2-PENTANOL | 6032-29-7 | HARMFUL/FLAMMABLE | SOLVENT | * | | |
| 2-PENTANONE | 107-87-9 | HARMFUL/FLAMMABLE | SOLVENT | * | | |
| 2-PROPANOL | 67-63-0 | FLAMMABLE | ALCOHOL | * | | |
| 3,4-DICHLOROBUTENE | | TOXIC/IRRITANT/FLAMMABLE | | * | | 96,04,15,001 |
| 30952 SYNTHOPON DAH94-505 | | IRRITANT (SEE STYRENE) | | * | | 01,01,23,006 |
| 3-AMINO-1-PROPANOL | 156-87-6 | SEE PROPANOLAMINE-1,3 | | | | |
| 3-AMINO-3,5,5-TRIMETHYLCYCLOHEXYLAMINE | 2855-13-2 | SEE ISOPHORENEDIAMINE | | | | |
| 3-AMINOPROPYLTRIETHOXYSIANE | 919-30-2 | CORROSIVE | BASE | * | | 98,09,09,001 |
| 3D TRASAR 3DT134 | | | WATER TREATMENT | * | | 070810004/MSDS |
| 3-DIETHYLAMINO-1-PROPYLAMINE | 104-78-9 | CORROSIVE / HARMFUL/FLAMMABLE | | * | | |
| 3-DIMETHYLAMINO-1-PROPYLAMINE | 109-55-7 | HARMFUL / CORROSIVE/FLAMMABLE | | * | | 95,05,29,013 |
| 3-HEPTANOL | 589-82-2 | HARMFUL/IRRITANT | SOLVENT | * | | |
| 3-METHYL-1-BUTANOL | 123-51-3 | HARMFUL/FLAMMABLE | SOLVENT | * | | |
| 3-OCTANOL | 589-98-0 | HARMFUL/IRRITANT | ALCOHOL | * | | |
| 3-PENTANOL | 584-02-1 | HARMFUL/FLAMMABLE | SOLVENT | * | | |
| 3-PENTANONE | 96-22-0 | FLAMMABLE | SOLVENT | * | | |
| 4,4-BIS-SULFOSTYRYL-BIPHENYL SODIUM SALT | 27344-41-8 | IRRITANT | | * | | 02,11,28,001/FDS |
| 4,4'-METHYLENEBIS(2-METHYLCYCLOHEXYLAMINE) | 6864-37-5 | TOXIC/CORROSIVE | | * | | |
| 4,4'-METHYLENEDICYCLOHEXYLDIISOCYANATE | 5124-30-1 | SEE DESMODUR W | | | | |
| 404a | | | HFC | * | | 99,06,14,006 |
| 407C | | | HFC | * | | 99,06,14,003 |
| 4-AMINOACETANILIDE | 122-80-5 | | | * | | 98,08,17,047 |
| 4-BROMOBENZYL BROMIDE | 589-15-1 | CORROSIVE | | * | | 03,08,06,001/FDS |
| 4-CHLORO-3-METHYLCRESOL | 59-50-7 | HARMFUL/DANGEROUS TO THE ENVIRONMENT | | * | | |
| 4-CHLORO-meta-CRESOL | 59-50-7 | SEE 4-CHLORO-3-METHYLCRESOL | | | | |
| 4-CHLOROPHENOL | 106-48-9 | HARMFUL/CORROSIVE | | * | | |
| 4-DIMETHYL AMINOBENZALDEHYDE | 100-10-7 | IRRITANT | | * | | 95,05,30,024 |
| 4-FLUORO-3-NITROANILINE | 364-76-1 | IRRITANT / HARMFUL | IRRITANT / NOXIOUS | * | | 95,06,01,017 |
| 4-FLUOROPHENOL | 371-41-5 | IRRITANT / HARMFUL | | * | | 99,11,18,001 |
| 4-METHOXY-2-NITROANILINE | 96-96-8 | TOXIC | | * | | 98,08,17,043 |
| 4-METHYL-2-PENTANONE | 108-10-1 | FLAMMABLE | | * | | 95,05,30,014 |
| 4-METHYLPENTAN-2-OL | 108-11-2 | IRRITANT,FLAMMABLE | SOLVENT | * | | 090330001/MSDS |
| 4-NITROBENZYL BROMIDE | 100-11-8 | CORROSIVE | | * | | 04,11,26,001/FDS |
| 4-OCTANOL | 589-62-8 | | SOLVENT | * | | |
| 5% LITHIUM TETRAHYDIDOALUMINATE IN THF | 16853-85-3 | CORROSIVE, Reacts with water | REDUCING AGENT | * | | 100624001/MSDS |
| 5% PALLADIUM/CHARBON CATALYST | | IRRITANT | | * | | 95,06,01,016 |
| 5-AMINO-ORTHO-CRESOL | 2835-95-2 | IRRITANT | | * | | 95,05,30,034 |
| 60% SODIUM HYDRIDE IN PARAFFINE OIL | 7646-69-7 | CORROSIVE, Reacts with water | BASE | * | | 100617001/MSDS |
| ABATE 100E INSECTICIDES | | TOXIC, IRRITANT | INSECTICIDE | * | | 100511003/MSDS |
| ABATE 105G MOSQUITO LARVICIDE GRANULES | | TOXIC, IRRITANT | INSECTICIDE SOLID | * | | 100511002/MSDS |
| ABATE 505G MOSQUITO LARVICIDE GRANULES | | TOXIC, IRRITANT | INSECTICIDE SOLID | * | | 100511001/MSDS |
| ABSOLU MOUSSE CHENE | | | | * | | 03,05,16,003/FDS |
| ABSOLU MOUSSE D'ARBRE | | | | * | | 03,05,16,002/FDS |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | | |
|--|------------|--|-----------------------|---|---|-------------------|
| ABSOLUE ALGUE MOUSSE ARBRE SUR HUILE VEGETALE E_3104 | 68648-41-9 | IRRITANT | | * | | 03,05,16,001/FDS |
| ACENAPHTHENE | 83-32-9 | IRRITANT | | * | | 95,06,19,004 |
| ACETALDEHYDE | 75-07-0 | HARMFUL/FLAMMABLE | ALDEHYDE | * | | 99,03,30,001 |
| ACETIC ACID | 64-19-7 | CORROSIVE/FLAMMABLE | ACID | * | | |
| ACETIC ACID POTASSIUM SALT | 127-08-2 | IRRITANT | ACID SALT | * | | 98,08,17,003 |
| ACETIC ACID SODIUM SALT | 6131-90-4 | IRRITANT | ACID SALT | * | | 98,08,17,004 |
| ACETIC ANHYDRIDE | 108-24-7 | CORROSIVE | | * | | |
| ACETONE | 67-64-1 | FLAMMABLE | SOLVENT | * | | |
| ACETONECYANHIDRINE | 75-86-5 | SEE 2-HYDROXY-2-METHYLPROPIONITRILE | | * | | 100414001/MSDS |
| ACETONITRILE | 75-05-8 | TOXIC FLAMMABLE | SOLVENT | * | | |
| ACETOPHENONE | 98-86-2 | HARMFUL | SOLVENT | * | | 99,07,02,004 |
| ACETYL CHLORIDE | 75-36-5 | CORROSIVE/FLAMMABLE/REACTS WITH WATER | ACID CHLORIDE | * | | 98,08,17,019 |
| ACETYLACETONE | 123-54-6 | HARMFUL / FLAMMABLE | | * | | 95,06,19,002 |
| ACID CHLORIDES | | | | * | | |
| ACID DICHLOROISOCYANURIC, SODIUM SALT | 2893-78-9 | SEE TROCLOSENE SODIUM | | * | | |
| ACID SOLUBLE OIL | | TOXIC, CORROSIVE | CONTENTS HF | | * | 071026001-MSDS |
| ACMOSAN 82-6001 | | IRRITANT | SOLVENT | * | | 01,01,23,004 |
| ACMOSOL 130-11 | | IRRITANT | SOLVENT | * | | |
| ACNOS 30-1003-111 | | IRRITANT | SOLVENT | * | | |
| ACQ 1900 | | HARMFUL/CORROSIVE/dangerous to the environment | BASE | * | | 070607001/FDS |
| ACROLEIN | 107-02-8 | SEE ACRYLALDEHYDE | ALKYLATING AGENT | * | | 071130001-MSDS |
| ACRYLALDEHYDE | 107-02-8 | VERY FLAMMABLE, VERY TOXIC, CORROSIVE | ALKYLATING AGENT | * | | 071130001-MSDS |
| ACRYLAMIDE | 79-06-1 | TOXIC | | * | | 95,06,19,007 |
| ACRYLIC ACID | 79-10-7 | CORROSIVE | ACID/ALKYLATING AGENT | * | | 98,08,17,006 |
| ACRYLONITRILE | 107-13-1 | TOXIC/FLAMMABLE | | * | | 070626001/MSDS |
| ACTICARBONEACIDE | | FLAMMABLE/IRRITANT | | * | | |
| ACTICATEUR N°2 | | CORROSIVE | ACID SALT | * | | |
| ACTICIDE 14 | | CORROSIVE | | * | | |
| ACTICIDE SPX | | IRRITANT | | * | | |
| ACTISIDE RS | | TOXIC, CORROSIVE, SENSITIZER | BIOCIDE | * | | 080225002/MSDS |
| ACTIVE COAL IN POWDER | | FLAMMABLE/IRRITANT | | * | | 95,06,01,019 |
| ADIPIC ACID | 124-04-9 | IRRITANT | CARBOXYLIC DIACID | * | | 98,08,17,007 |
| AGATE VILLE-ROUTE | | FLAMMABLE, IRRITANT, SENSITIZER | ACRYLATES | * | | 080326005/FDS |
| AGI GAZON EV | | HARMFUL | | * | | 05,01,20,016/FDS |
| AIBN | 78-67-1 | HARMFUL/EXPLOSIVE | | * | | 95,06,01,023 |
| AL ELECTROLYTE | | TOXIC/CORROSIVE | | * | | 060911001/MSDS |
| ALCOHOLIC POTASH | | CORROSIVE / TOXIC | | * | | 95,06,26,001 |
| ALD C-8 | 124-13-0 | IRRITANT | | * | | 03,05,16,011/FDS |
| ALD CYCLAMEN | 103-95-7 | | | * | | 03,05,16,007/FDS |
| ALDEHYDE C9 NONYLIC FOOD GRADE | 124-19-6 | IRRITANT | | * | | 03,05,16,006/FDS |
| ALGENOL SP | | IRRITANT | | * | | 03,05,16,008/FDS |
| ALGUE MOUSSE D'ALGUE CBA 9810/A | | IRRITANT | | * | | 03,05,16,009/FDS |
| ALKALI | | CORROSIVE | STRONG BASE | * | | |
| ALLYL ALCOHOL | 107-18-6 | TOXIC | SOLVENT | * | | |
| ALLYL AMYLGLYCOLATE | | HARMFUL | | * | | 03,05,16,010/FDS |
| ALOBIO 492 | | | | * | | 01,05,20,002/FDS |
| ALODINE 1200 | | TOXIC/CORROSIVE | | * | * | 03,09,18,002/FDS |
| ALODINE 4780 | | TOXIC/CORROSIVE | | * | * | |
| ALOSUC 845D | | HARMFUL | | * | | 05,01,20,017/FDS |
| ALOSUC 854A | | | | * | | 01,05,20,003/FDS |
| ALPHA-D-LACTOSE | 14641-93-1 | | | * | | 04,12,07,006/MSDS |
| alpha-METHYLSTYRENE | 98-83-9 | IRRITANT | AROMATIC | * | | 99,07,02,003 |
| ALPHASET RESIN TPA85 | | SENSITIZER | | * | | 04,04,06,004/FDS |
| ALS | 10043-01-3 | IRRITANT | | * | | 02,02,19,002 |
| ALUMINIUM FLUORIDE | 7784-18-1 | HARMFUL | SOLID | * | | 071127001/NEP |
| ALUMINIUM HYDROXIDE CHLORIDE | 1327-41-9 | CORROSIVE | | * | | 02,02,19,003 |
| ALUMINIUM SULFATE | 10043-01-3 | IRRITANT | | * | | 02,02,19,002 |
| ALUMINIUM SULFATE Hexadecahydrate | 16828-11-8 | | ALUMINIUM SALT | * | | |
| AMBERLITE IR-120 Resin | | | | * | | 95,05,30,016 |
| AMBERLITE IRN 150 L | | IRRITANT | | * | | 95,05,30,038 |
| AMIDES | | | | * | | |
| AMINES | | | | * | | |
| AMINO-2 ETHANOL | 141-43-5 | HARMFUL/IRRITANT | ORGANIC BASE | * | | |
| AMINO-4 NITRO-3-PHENOL | 610-81-1 | HARMFUL/IRRITANT | | * | | 95,05,30,005 |
| AMINOACETIC ACID | 56-40-6 | SEE GLYCINE | | * | | |
| AMMEDIOL | 115-69-5 | SEE 2-AMINO-2-METHYL-1,3-PROPANEDIOL | | * | | |
| AMMONIA | 7664-41-7 | IRRITANT GAZ | GAS | * | | |
| AMMONIUM BICARBONATE | 1066-33-7 | HARMFUL | AMMONIA SALT | * | | |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | |
|--|------------|--|-------------------|---|-------------------|
| AMMONIUM BIFLUORIDE | 1341-49-7 | TOXIC/CORROSIVE | FLUORIDE | * | |
| AMMONIUM CHLORIDE | 12125-02-9 | HARMFUL | AMMONIA SALT | * | |
| AMMONIUM FLUORIDE 40% | 12125-01-8 | TOXIC | | * | 06.03.02.001/MSDS |
| AMMONIUM FLUORIDE ETCHANT AF 87,5-12,5 VLSI SELECTIPUR | | | | | 04,06,24,002/MSDS |
| AMMONIUM HEPTAMOLYBDATE | 12054-85-2 | HARMFUL/IRRITANT | | * | |
| AMMONIUM HYDROGEN FLUORIDE | 1341-49-7 | TOXIC/CORROSIVE | FLUORIDE | * | 04,03,17,001/PDT |
| AMMONIUM HYDROXIDE | 1336-21-6 | CORROSIVE/DANGEROUS TO THE ENVIRONMENT | BASE | * | |
| AMMONIUM MOLYBDATE TETRAHYDRATE | 12054-85-2 | SEE AMMONIUM HEPTAMOLYBDATE | | | |
| AMMONIUM NITRATE | 6484-52-2 | IRRITANT/OXIDIZING | AMMONIA SALT | * | |
| AMMONIUM SULFATE | 7783-20-2 | | AMMONIA SALT | * | |
| AMMONIUM THIOCYANATE | 1762-95-4 | HARMFUL/IRRITANT | | * | 01,12,04,004 |
| AMP | 124-68-5 | SEE 2-AMINO-2-METHYL-1-PROPANOL | | | |
| AMPD | 115-69-5 | SEE 2-AMINO-2-METHYL-1,3-PROPANEDIOL | | | |
| AMYL ACETATE | 628-63-7 | FLAMMABLE | ESTER | * | |
| ANHYDROUS ALUMINIUM CHLORIDE | 7446-70-0 | CORROSIVE | ACID SALT | * | |
| ANHYDROUS COBALT CHLORIDE | 7646-79-9 | TOXIC | | * | 95,06,26,008 |
| ANILINE | 62-53-3 | TOXIC/DANGEROUS TO THE ENVIRONMENT | BASE | * | |
| ANTHRACENIC OIL | | HARMFUL | NOXIOUS | * | |
| ANTOX 71 E | | CORROSIVE/TOXIC | MIX HF/HNO3 | * | 99,09,13,001 |
| AQUAREX 2290 | | HARMFUL | | * | 03,11,24,003/FDS |
| ARALDITE | | | RESIN | * | |
| ARCIE 2111 | | CORROSIVE | ZINC SALT | * | 03,11,24,004/FDS |
| ARDROX 2106R | | HARMFUL/IRRITANT | SOLVENT | * | 070326001/MSDS |
| ARDROX 2480 | | HARMFUL/CORROSIVE | | * | 02,07,11,002 |
| ARDROX LEEDER 1900B | | IRRITANT | BASE | * | 99,09,07,001 |
| ARMEEN 2HD | 61789-79-5 | IRRITANT | BASE | * | 01,01,10,004 |
| ARMEEN 2HT | 61789-79-5 | HARMFUL/IRRITANT | AMINE | * | 061218001/EXP |
| ARMEEN 2MHTD | 61788-95-2 | CORROSIVE/DANGEROUS TO THE ENVIRONMENT | BASE | * | 01,01,10,002 |
| ARMEEN 2MMCD | | CORROSIVE | BASE | * | 01,01,10,005 |
| ARMEEN C | 61788-46-3 | CORROSIVE | BASE | * | 01,01,10,001 |
| ARMEEN CD | 61788-46-3 | HARMFUL/CORROSIVE | AMINE | * | 061218001/EXP |
| ARMEEN M | 68037-92-3 | IRRITANT/DANGEROUS TO THE ENVIRONMENT | AMINES | * | |
| ARMEEN M2HT | 61788-63-4 | HARMFUL/IRRITANT | AMINE | * | 061218001/EXP |
| ARMEEN TD | 61790-33-8 | HARMFUL/CORROSIVE | AMINE | * | 061218001/EXP |
| AROPOL TM Di4s253 | | IRRITANT (SEE STYRENE) | | * | 01,01,23,007 |
| ARQUAD MCB-50 | | HARMFUL | | * | 01,01,10,008 |
| ARSENIC ACID | 7778-39-4 | TOXIC/CORROSIVE | ACID | * | 99,02,03,001 |
| ARYLSULFONIC ACIDS | | CORROSIVE | | * | |
| AS 140 | | TOXIC | | * | |
| ASCORBIC ACID | 50-81-7 | | WEAK ACID | * | |
| ASPHALT | | HARMFUL/FLAMMABLE | | * | |
| AUROCA RENOVATEUR INOX | | | | * | 060504001/FDS |
| AXION PS 3430 | | CMR, TOXIC, IRRITANT | ALKYLATING AGENT | * | 090506001/MSDS |
| AZ 100 REMOVER | | | | * | 04,06,24,001/MSDS |
| AZ 400 K DEVELOPER | | | | * | 04,06,24,001/MSDS |
| AZOBISOBUTYRONITRILE | 78-67-1 | SEE AIBN | | * | |
| BACTOLYSE 24 | | | | * | 89,09,04,006 |
| BACTOLYSE 66 | | | | * | 89,09,04,007 |
| BAKER REZITM-28 DILUTE | | | | * | 04,06,24,001/MSDS |
| BAKER REZITM-38 CONCENTRATED | | | | * | 04,06,24,001/MSDS |
| BANDAX ENDUIT | | IRRITANT, SENSITIZER | ACRYLATES | * | 080326004/FDS |
| BARDHAL | | IRRITANT | | * | |
| BARIUM CHLORIDE | 10326-27-9 | HARMFUL | | * | 95,06,26,007 |
| BASES | | | | * | |
| BAUXITE SLURRY | | CORROSIVE | BASE | * | 01,10,15,001 |
| BAYER LIQUOR | | CORROSIVE | BASE | * | 99,02,05,001 |
| BENZALDEHYDE | 100-52-7 | HARMFUL | AROMATIC ALDEHYDE | * | 99,08,17,017 |
| BENZALKONIUM CHLORIDE | 8001-54-5 | CORROSIVE | | * | 04,04,13,001/FDS |
| BENZENE | 71-43-2 | FLAMMABLE TOXIC | SOLVENT | * | |
| BENZENE-1,4-DICARBOXALDEHYDE | 623-27-8 | SEE TEREPHTHALALDEHYDE | | | |
| BENZOIC ACID | 65-85-0 | HARMFUL | WEAK ACID | * | |
| BENZONITRIL SYNTHESSE KWALITEIT | 100-47-0 | IRRITANT (SEE BENZONITRILE) | | * | 01,12,04,002 |
| BENZONITRILE | 100-47-0 | IRRITANT | | * | 01,12,04,002 |
| BENZOYL PEROXIDE 75 % IN WATER | 94-36-0 | TOXIC/IRRITANT | | * | 95,06,01,020 |
| BENZYL ALCOHOL | 100-51-6 | HARMFUL SOLVENT | SOLVENT | * | |
| BENZYL CHLORIDE | 100-44-7 | IRRITANT | | * | |
| BENZYLAMINE | 100-46-9 | HARMFUL/CORROSIVE | | * | 02,07,05,003 |
| BENZYLcyanide | 140-29-4 | HARMFUL/IRRITANT | ACID | * | 01,12,04,003 |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | | |
|--|------------|---|---------------------|---|---|-------------------|
| BENZYL CYANIDE SYNTHESIS KVALITEIT | 140-29-4 | HARMFUL/IRRITANT (SEE BENZYL CYANIDE) | ACID | * | | 01,12,04,003 |
| BENZYL DIMETHYLAMMONIUM CHLORIDE | 61789-71-7 | HARMFUL | | * | | 01,01,10,008 |
| BEVEDAN 1 | | SEE DIPHENYLMETHANE DIISOCYANATE PREPOLYMERISED | PREPOLYMER | * | | 100315002/MSDS |
| BEVEDOL S21 | | | POLYURETHAN RESIN | * | | 100315001/MSDS |
| BEVEDOL WF21 | | | POLYURETHAN RESIN | * | | 100315001/MSDS |
| BHF 39/1 | | | | | * | 06,03,02,003/MSDS |
| BIOCIDE 195 EX | | TOXIC/CORROSIVE | | * | | 03,11,24,002/FDS |
| BIO-KAT-1006-WU | | | | * | | 04,07,01,001 |
| BIOSPERSE 250 | | CORROSIVE, SENSITIZER | | * | | 071122001-MSDS |
| BIS-(2-HYDROXYPROPYL) AMINE | 110-97-4 | SEE DIISOPROPANOLAMINE | | | | |
| BIS(a,a-DIMETHYLBENZYL) PEROXIDE | 80-43-3 | OXYDISING, IRRITANT | OXYDIZER | * | | 100413001/MSDS |
| BISTRIFLUOROACETAMIDE | 407-24-9 | IRRITANT / FLAMMABLE | | * | | 95,06,26,009 |
| BITUME | | HARMFUL/FLAMMABLE | | * | | |
| BLASOCUT KOMBI | | | EMULSIVE LIQUID | * | | 06,06,14,002/MSDS |
| BLEACH | 7681-52-9 | CORROSIVE | BASE | * | | |
| BORANE TETRAHYDROFURANE COMPLEX SOLUTION | | CORROSIVE, Reacts with water | REDUCING AGENT | * | | 100617002/MSDS |
| BORAX | 1303-96-4 | SEE SODIUM TETRABORATE DECAHYDRATE | | | | |
| BORDEN RESIN UL200 | | | | * | | 04,12,15,001/FDS |
| BORIC ACID | 10043-35-3 | HARMFUL | ACID | * | | |
| BOROL | | HARMFUL/CORROSIVE | BASE/REDUCING AGENT | * | | 02,03,25,001 |
| BORON NITRIDE SOLID | | | | * | | 04,04,06,015/MSDS |
| BORON TRICHLORIDE | 10294-34-5 | TOXIC, Reacts with water | ACID | * | | 100816001/MSDS |
| BORON TRIFLUORIDE ANHYDROUS | 2095581 | TOXIC/CORROSIVE/REACTS WITH WATER | | | * | |
| BORON TRIFLUORIDE DIHYDRATE | 13319-75-0 | TOXIC/CORROSIVE | | | * | |
| BORON TRIFLUORIDE DIMETHYL ETHERATE | 353-42-4 | CORROSIVE/FLAMMABLE/REACTS WITH WATER | | | * | |
| BORON TRIFLUORIDE ETHERATE | 109-63-7 | CORROSIVE/FLAMMABLE/REACTS WITH WATER | | | * | |
| BORON TRIFLUORIDE IN ACETIC ACID | 373-61-5 | CORROSIVE/HARMFUL | | | * | 95,09,12,001 |
| BOUZE0 MB 271 | | | | * | | 06,03,01,003/FDS |
| BRAKE OIL POCLAIN | | | | * | | |
| BROMOACETIC ACID | 79-08-3 | TOXIC/CORROSIVE | ACID | * | | 98,08,17,008 |
| BROMOCRESOL PURPLE | 115-40-2 | | | * | | 95,06,21,004 |
| BROMODODECANE | 143-15-7 | | | * | | |
| BROMOFORM | 75-75-2 | TOXIC/IRRITANT/DANGEROUS TO THE ENVIRONMENT | | * | | |
| BROMOMETHANE | 74,83,9 | TOXIC | | * | | 03,04,01,002/FDS |
| BROMOTRIMETHYLSILANE | 2857-97-8 | CORROSIVE | | * | | 00,10,23,001 |
| BTFA | 407-24-9 | SEE BISTRIFLUOROACETAMIDE | | | | |
| BUFFALO SH | | | | * | | 05,01,20,011/FDS |
| BUFFER SOLUTIONS 4,7,9,10 | | | | * | | 95,06,26,005 |
| BULCOAT 33 | 17439-11-1 | TOXIC/CORROSIVE | | | * | 05,04,05,001/FDS |
| BUSAN 881 | | HARMFUL, CORROSIVE, SENSITIZER | BIOCID | * | | 080222003/MSDS |
| BUTANAL | 123-72-8 | CORROSIVE/SENSITIZER | REDUCING AGENT | * | | 070927004/MSDS |
| BUTANOL-1 | 71-36-3 | HARMFUL/FLAMMABLE | SOLVENT | * | | |
| BUTANOL-2 | 14898-79-4 | HARMFUL/FLAMMABLE | SOLVENT | * | | |
| BUTANOX | | | OXIDIZER | * | | |
| BUTENONE | 78-94-4 | TOXIC/CORROSIVE | | * | | 070927001/MSDS |
| BUTYL ACETATE | 123-86-4 | FLAMMABLE | SOLVENT | * | | |
| BUTYL NITRITE | 544-16-1 | TOXIC/FLAMMABLE | | * | | 95,06,01,018 |
| BUTYLACRYLATE | 141-32-2 | IRRITANT/SENSITIZER | | * | | 070927002/MSDS |
| BUTYLDIGLYCOL | 112-34-5 | SEE DIETHYLENE GLYCOL MONOBUTYLEETHER | | | | |
| BUTYLENEGLYCOL | | SEE BUTANE DIOL | | | | |
| BUTYLGLYCOL | 111-76-2 | SEE 2-BUTOXYETHANOL | | | | |
| BUTYLGLYCOL ACETATE | 112-07-2 | HARMFUL | | * | | |
| BUTYRALDEHYDE | 123-72-8 | SEE BUTANAL | | * | | |
| BUTYROLACTONE | 96-48-0 | HARMFUL/IRRITANT | SOLVENT | * | | 070716001/MSDS |
| CAB-O-SIL Untreated Fumed Silica | | | | | | 04,12,07,009/FDS |
| CaD 11,21,23 | | HARMFUL, IRRITANT | | * | | 080909001/MSDS |
| CADMIUM SULPHIDE | 1306-23-6 | CMR, TOXIC | | * | | 100927001/MSDS |
| CAFFEINE | 58-08-2 | HARMFUL | | * | | |
| CALCIUM CARBONATE | 471-34-1 | | CALCIUM SALT | * | | |
| CALCIUM CHLORIDE (anhydrous) | 10043-52-4 | IRRITANT | CALCIUM SALT | * | | |
| CALCIUM HYDROXIDE | 1305-62-0 | IRRITANT/CORROSIVE | IRRITANT | * | | |
| CALCIUM HYPOCHLORITE | 7778-54-3 | CORROSIVE/OXIDIZING | OXIDIZER | * | | |
| CALCIUM NITRATE TETRAHYDRATE | 13477-34-4 | IRRITANT/OXIDIZING | | * | | |
| CALCIUM OXIDE | 1305-78-8 | IRRITANT | BASE | * | | |
| CALOXOL W3 | | IRRITANT | DEHYDRATING | * | | |
| CAMPHOR | 21368-68-3 | IRRITANT/NOCIF/FLAMMABLE | | * | | 95,05,30,012 |
| CAMPHOSULFONIC ACID RACEMIC | 1450959 | CORROSIVE | | * | | 95,05,30,028 |
| CANYON | | | | * | | 05,01,20,012/FDS |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | | |
|---------------------------------------|------------|--|--------------------|---|---|-------------------|
| CAPRIC ACID | 334-48-5 | IRRITANT | | * | | 95,05,30,008 |
| CAPROLACTAM | 105-60-2 | HARMFUL/IRRITANT | SOLVENT | * | | 070716002/MSDS |
| CARBARYL | 63-25-2 | HARMFUL/IRRITANT | SOLID | * | | 070301001/MSDS |
| CARBIDE CALCIUM | 75-20-7 | REACTS WITH WATER | | | | |
| CARBOLIC OIL | 84650-03-3 | TOXIC/CORROSIVE | | * | | 04,11,26,003/FDS |
| CARBON DISULFIDE | 75-15-0 | TOXIC/FLAMMABLE | SOLVENT | * | | |
| CARBON TETRACHLORIDE | 56-23-5 | VERY TOXIC/DANGEROUS TO THE ENVIRONMENT | | * | | 95,06,26,003 |
| CARBOSULF | | CORROSIVE (SEE SODIUMHYDROSULFITE) | | * | | 01,02,08,001 |
| CARBOTHIX 1:2 KOMP.A | | IRRITANT | BASE | * | | 100316002/MSDS |
| CARBOTHIX 1:2 KOMP.B | | HARMFUL, IRRITANT, SENSITIZER | | * | | 100316001/MSDS |
| CARBUREACTOR JET A1 | 8008-20-6 | FLAMMABLE | HYDROCARBURES MIX | * | | 99,02,22,002 |
| CATALYSEUR 3052 JEANNEAU | | | | * | | |
| CAUSTIC SODA LIQUOR | | CORROSIVE | STRONG BASE | * | | 070810005/MSDS |
| CDR 1045 | | CORROSIVE/TOXIC | MIX HF/HNO3 | | * | 99,02,01,001 |
| CELCURE CCA TYPE C50 | | | | * | | 06.04.10.003/FDS |
| CELECURE AC - 800 | | | | * | | 06.04.10.004/FDS |
| CELLOBOND J2027L | | HARMFUL/TOXIC | | * | | 04,04,06,005/FDS |
| CELLOSIZ | | SEE HYDROXYETHYLCELLULOSE | | | | |
| CERAMEX | | | PICKLING AGENT | | * | 05,06,23,002/MSDS |
| CERTINCOAT TC 100 | 1118-46-3 | CORROSIVE | | * | | 02,08,30,003 |
| CETANE IMPROVER | | SEE 2-ETHYLHEXYL NITRATE | | * | | |
| CETYL ALCOHOL | 36653-82-4 | SEE 1-HEXADECANOL | | | | |
| CHALOXYD | | | STRONG OXIDIZER | * | | |
| CHEMIUM ETCH LSI SELECTIPUR | | | | * | | 04,06,24,001/MSDS |
| CHITOSAN | 9012-76-4 | | | * | | 95,05,30,025 |
| CHLORINE | 7782-50-5 | TOXIC/IRRITANT/DANGEROUS TO THE ENVIRONMENT | GAS | * | | |
| CHLORINE DIOXIDE 2-2,7% | | VERY TOXIC, IRRITANT, Dangerous to the environment | STRONG OXYDISER | * | | 080618001/MSDS |
| CHLORO ISOTHIAZOLINE | | IRRITANT | | * | | |
| CHLOROACETAMIDE | 79-07-2 | TOXIC / IRRITANT | | * | | 95,05,29,011 |
| CHLOROACETIC ACID | 79-11-8 | TOXIC/DANGEROUS TO THE ENVIRONMENT | | * | | 070712001/MSDS |
| CHLOROACETYL CHLORIDE | 79-04-9 | TOXIC/CORROSIVE/REACTS WITH WATER | | * | | 06,03,01,002/MSDS |
| CHLOROBENZENE | 108-90-7 | HARMFUL/FLAMMABLE/DANGEROUS TO THE ENVIRONMENT | | * | | 95,05,30,017 |
| CHLOROFORM | 67-66-3 | HARMFUL/IRRITANT | | * | | 95,06,19,012 |
| CHLOROMETHYL-ISOTHIAZOLINONE | 26172-55-4 | CORROSIVE | | * | | 01,12,04,001 |
| CHLOROSULFONIC ACID | 7790-94-5 | CORROSIVE | STRONG ACID | * | | 98,08,17,009 |
| CHLOROTRIMETHYLSILANE | 75-77-4 | SEE TRIMETHYLCHLOROSILANE | | | | |
| CHOLESTEROL | 57-88-5 | | | * | | 95,05,29,001 |
| CHROMIC ACID | 1333-82-0 | TOXIC/CORROSIVE/OXIDIZING | | * | | |
| CINNAMYL ALCOHOL | 104-54-1 | HARMFUL/IRRITANT | ALCOHOL | * | | |
| CIS-1,4-DICHLORO-2-BUTENE | 1476-11-5 | TOXIC/IRRITANT/FLAMMABLE | | * | | 96,04,15,001 |
| CIS-9-OCTADECEN-1-OL | 143-28-2 | SEE OLEIC ALCOHOL | | | | |
| CIS-9-OCTADECENOIC ACID | 112-80-1 | SEE OLEIC ACID | | | | |
| CIT 9387 | | | DYEING IN SOLUTION | * | | |
| CITGO GASOLINES | | | HYDROCARBURES MIX | * | | 061218001/MSDS |
| CITRIC ACID | 77-92-9 | IRRITANT | WEAK ACID | * | | |
| CITROFOL B II | 77-90-7 | | SOLVENT | * | | 05,01,20,016/MSDS |
| CLEANSEPT UBI | | CORROSIVE, SENSITIZER | BIOCIDE | * | | 080222002/MSDS |
| CN, CR and CS SPRAY or GAS | | IRRITANT | | * | 05,08,30,001/biblio-tests and publication Hall AH J Emerg M | |
| COBALT CHLORIDE HEXAHYDRATE | 7791-13-1 | HARMFUL | | * | | 95,06,26,008 |
| COCOALKYL AMINE | | CORROSIVE | | * | | 95,06,01,002 |
| COMORAL DP | | IRRITANT | SKIMMING PRODUCT | * | | 99,11,02,005 |
| concentrated n-BUTYLLITHIUM in HEXANE | | SEE n-BUTYLLITHIUM | | * | | 071004004-MSDS |
| COPPER SULFATE ANHYDROUS | 7758-98-7 | HARMFUL/IRRITANT | | * | | |
| COPRAH FAT ACID | | IRRITANT | | * | | 95,06,01,012 |
| CORTROL OS5601 | | IRRITANT, SENSITIZER | REDUCING AGENT | * | | 100310001/MSDS |
| CPR 154 F (PAAE Resin) | | see polyamidoamine-epichlorhydrine (PAAE) resin in water | | * | | 080225003/MSDS |
| Créosote OIL WEI type C | | | | * | | 06.04.10.002/FDS |
| CROTONIC ACID | 107-93-7 | HARMFUL / CORROSIVE | | * | | 95,05,30,018 |
| CRYOLITHE | 15096-52-3 | TOXIC/CORROSIVE AT MELTING POINT | | * | | |
| CRYSOAT® OC | | CMR,HARMFUL, CORROSIVE | ACID OXYDIZER | * | | 100525001/MSDS |
| CUMENE | 98-82-8 | IRRITANT | AROMATIC | * | | 99,07,02,001 |
| CUMENE HYDROPEROXYDE | 80-15-9 | FLAMMABLE CORROSIVE/OXIDIZING | PEROXIDE | * | | 91,03,06,001 |
| CUPOSIT ACCELERATOR 19H | | IRRITANT | | * | | 061026001/MSDS |
| CYANOACRYLATE | | IRRITANT/TOXIC | | * | | 96,05,06,001 |
| CYCLOHEXANE | 110-82-7 | FLAMMABLE | | * | | 95,06,19,014 |
| CYCLOHEXANONE | 108-94-1 | HARMFUL/FLAMMABLE | | * | | |
| CYCLOHEXYLAMINE | 108-91-8 | CORROSIVE/HARMFUL/FLAMMABLE | | * | | 95,02,09,001 |
| CYSTEIC ACID | 23537-25-9 | CORROSIVE | ACID | * | | 98,08,17,010 |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | | |
|--|------------|---------------------------------------|-------------------|---|---|-------------------------|
| DALTOFLEX 20650 | | | | * | | 95,03,15,003 |
| DALTOFLEX X ES 11000 | | | | * | | 95,03,15,001 |
| DARTOLINE SRB 71 | | | SKIMMING PRODUCT | * | | 89,09,04,002 |
| DBTCI | 683-18-1 | SEE DIBUTYLTIN DICHLORIDE | | * | | |
| DCC | 538-75-0 | TOXIC | | * | | 03,02,25,003/FDS |
| DECAPANT CASTOLIN | | HARMFUL | | * | | 04,02,19,001/MSDS |
| DECYL ALCOHOL | 112-30-1 | SEE DECANOL | | * | | |
| DEF-AD-80 WU | | | | * | | 04,07,01,003/MSDS |
| DEGAFAS | 1758-73-2 | IRRITANT | | * | | 02,03,04,001 |
| DEGREASER O3 | | IRRITANT | ALKALIN CLEANSER | * | | 080618002/MSDS |
| DELAN WG | | SEE DITHIANON | FUNGICIDE | * | | 080728001/MSDS |
| DELTA-T-2593 | | | | * | | 90,01,22,015 |
| DELTA-T-2834 | | | | * | | 90,01,22,014 |
| DEOXIDINE SC 56 CF | | TOXIC/CORROSIVE | | * | * | |
| DESMODUR N-3200 | 28182-81-2 | TOXIC | POLYDIISOCYANATE | * | | 98,02,26,003 |
| DESMODUR W | 5124-30-1 | TOXIC | DIISOCYANATE | * | | 98,02,26,002 |
| DESMOPHEN 550 U | | | POLYOL POLYETHER | * | | 98,02,26,001 |
| DESMORAPID 726 B | | HARMFUL/CORROSIVE | | * | | 05,01,20,014/FDS |
| DEVELOPER E3/51 | | | | * | | 04,06,24,001/MSDS |
| DEX LO CL | | HARMFUL | | * | | |
| DIAMINE LAURYLPROPYLENE | 5538-95-4 | CORROSIVE | BASE | * | | 02,03,13,001 |
| DIBENZYLTOLEUENE | 26898-17-9 | | | * | | 02,12,10,001/FDS |
| DIBUTYLBIS(PENTANE-2,4-DIONATO-O,O')TIN | 22673-19-4 | IRRITANT | CATALYSEUR | * | | 080509002/MSDS |
| DIBUTYLTIN DICHLORIDE | 683-18-1 | CMR, TOXIC, CORROSIVE | ALKYLATING AGENT | * | | 090430001/MSDS |
| DICHLORO-1,1-FLUORO-ETHANE | 1717-00-6 | | | * | | 02,08,01,002 |
| DICHLOROMETHANE | 75-09-2 | HARMFUL | SOLVENT | * | | |
| DICHLOROSILANE | 4109-96-0 | FLAMMABLE/CORROSIVE/ TOXIC | GAS | * | | 04,04,06,013/MSDS |
| DICUMYLE PEROXIDE | 80-43-3 | SEE BIS(a, a-DIMETHYLBENZYL) PEROXIDE | | * | | 100413001/MSDS |
| DIESEL FUEL ADDITIVE AC1000 | | CMR, HARMFUL, IRRITANT, SENSITIZER | SOLVENT | * | | 100225002/MSDS |
| DIESEL FUEL ADDITIVE EMDFA 301 CDG | | CMR, HARMFUL, IRRITANT | SOLVENT | * | | 100225001/MSDS |
| DIESEL OIL | 68334-30-5 | HARMFUL/FLAMMABLE | HYDROCARBURES MIX | * | | 99,02,22,004 |
| DIETHANOLAMINE | 111-42-2 | IRRITANT | AMINE | * | | |
| DIETHYL ETHER | 60-29-7 | FLAMMABLE | | * | | |
| DIETHYL SULFATE | 64-67-5 | TOXIC | | * | | 04,02,05,003/FDS |
| DIETHYLAMINE | 109-89-7 | CORROSIVE/FLAMMABLE | BASE | * | | 98,08,17,028 |
| DIETHYLENE DIAMINE | 110-85-0 | SEE PIPERAZINE ANHYDROUS | | * | | |
| DIETHYLENE GLYCOL MONOBUTYL ETHER | 112-34-5 | IRRITANT | | * | | |
| DIETHYLENGLYCOL | 111-46-6 | HARMFUL | DIOL | * | | 95,05,29,010 |
| DIETHYLENETRIAMINE | 111-40-0 | CORROSIVE | BASE | * | | 98,08,17,029 |
| DIGRAIN CONCENTRE | | | INSECTICIDE | * | | 90,01,22,017 |
| DIISOPROPANOLAMINE | 110-97-4 | IRRITANT | | * | | 95,05,29,017 |
| DIISOPROPYL ETHER | 108-20-3 | FLAMMABLE | | * | | 95,06,19,010 |
| DIMETHYL C12/C14 FAT AMINE | | SEE DODECYLDIMETHYLAMINE | | * | | |
| DIMETHYL SULFATE | 77-78-1 | TOXIC / CORROSIVE / CARCINOGENIC | | * | | 95,05,30,015 |
| DIMETHYL SULPHIDE | 75-18-3 | IRRITANT/FLAMMABLE | | * | | 070903001/MSDS |
| DIMETHYLAMINE, 40wt.% solution in water | 124-40-3 | HIGHLY FLAMMABLE/CORROSIVE | | * | | |
| DIMETHYLAMINOETHYL METHACRYLATE | 2867-47-2 | HARMFUL/IRRITANT | | * | | 95,06,01,015 |
| DIMETHYLENE OXIDE | 75-21-8 | SEE ETHYLENE OXIDE | | * | | |
| DIMETHYLFORMAMIDE | 68-12-2 | TOXIC/HARMFUL/IRRITANT | NOXIOUS | * | | |
| DIMETHYLMETHYLENEBIS(CYCLOHEXYLAMINE) | 6864-37-5 | TOXIC/CORROSIVE | | * | | 97,06,04,005 |
| DIMETHYLPHTALAT | 131-11-3 | IRRITANT | | * | | 01,01,23,008 |
| DIMETHYLSULFOXIDE | 67-68-5 | IRRITANT | | * | | 95,06,20,001 |
| DI-n-BUTYLPHTHALATE | 84-74-2 | IRRITANT | | * | | 03,06,04,001/MSDS |
| DINORAM C | | CORROSIVE | | * | | |
| DIPHENYLMETHANEDIISOCYANATE PREPOLYMERISED | 9016-87-9 | HARMFUL, IRRITANT, SENSITIZER | PREPOLYMER | * | | 100315003/MSDS |
| DIPHYL | | IRRITANT | | * | | 02,05,15,001/F-PROMETRA |
| DIPOTASSIUM HYDROGEN PHOSPHATE ANHYDROUS | 2139900 | | | * | | |
| DISODIUM HEXAFLUOROSILICATE | 16893-85-9 | TOXIC | ACID | * | * | 100505001/MSDS |
| DISODIUM HYDROGEN PHOSPHATE ANHYDROUS | 7558-79-4 | | | * | | |
| DISODIUM HYDROGEN PHOSPHATE DIHYDRATE | 10039-32-4 | | | * | | |
| DISODIUM HYDROGEN PHOSPHATE DIHYDRATE | 10028-24-7 | | | * | | |
| DISODIUM HYDROGEN PHOSPHATEHEPTAHYDRATE | 7782-85-6 | | | * | | |
| DISSOLVINE E-39-S | | CORROSIVE | BASE | * | | 01,01,10,006 |
| DITHIANON | 3347-22-6 | HARMFUL, IRRITANT | FUNGICIDE | * | | 080728001/MSDS |
| DL-SERINE | 302-84-1 | | AMINED ACID | * | | 98,08,17,032 |
| D-MANNITOL | 69-65-8 | | ALCOHOL | * | | |
| DMF | 68-12-2 | SEE DIMETHYLFORMAMIDE | | * | | |
| DMSO | 67-68-5 | SEE DIMETHYLSULFOXIDE | | * | | |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | | |
|---|------------|--|---|-----------------------|---|---------------------|
| DODECANAL | 112-54-9 | | | * | | 03,05,16,005/FDS |
| DODECYLBENZENE SULFONIC ACID | 27176-87-0 | | CORROSIVE | ACID | * | 02,02,19,005 |
| DODIFLOW 4134/ROAD TANKER | | | HARMFUL | SOLVENT AND POLYMERES | * | 090204003/MSDS |
| DODISTAT 5549 | | | CMR, HARMFUL, IRRITANT | SOLVENT AND POLYMERES | * | 090204009/MSDS |
| DOMESTOS | | | IRRITANT | BASE, OXIDIZER | * | 080931001/MSDS |
| DOMOPOL 5110 | | | IRRITANT | POLYESTER RESIN | * | 100823001/MSDS |
| DOW CORNING 561 SILICONE | | | | COOLING LIQUID | * | |
| DOWANOL 33-B | | | | | * | |
| DREWCLEAN C 506 | | | IRRITANT | | * | 04,04,06,001/MSDS |
| DREWCLEAN C 705 | | | IRRITANT | | * | 04,04,06,002/MSDS |
| D-SODIUM GLUCONATE | 527-07-1 | | | SODIUM SALT | * | |
| D-SORBITOL | 50-70-4 | | | SWEETENED | * | |
| DUCK BILE | | | CORROSIVE | CORROSIVE | * | 05,07,06,001 |
| DUPONT VELPAR® | | | HARMFUL, IRRITANT | PESTICIDE | * | 100512001/MSDS |
| DURCISSEUR 2564 | | | | | * | 04,12,15,002/FDS |
| DYNASOLVE 210 | | | CORROSIVE/CARCINOGENIC | SOLVENT | * | 070206001/MSDS |
| ECTIBAN EC | | | SENSITIZER/DANGEROUS TO THE ENVIRONMENT | SOLVENT | * | 070226001/MSDS |
| EDTA | 60-00-4 | | SEE ETHYLENE DIAMINE TETRAACETIQUE ACID | | | |
| ELBEFLEX RC 1 RH 648 A | | | IRRITANT | PAINT | * | 080326003/FDS |
| ENERPAC HF HYDRAULIC OIL | | | IRRITANT | VISCOUS OIL | * | 1001129001/MSDS |
| EPICHLORHYDRINE | 106-89-8 | | CARCINOGENIC/VERY TOXIC/CORROSIVE/FLAMMABLE | | * | 95,06,02,005 |
| EPOXY-1,2-HEXADECANE | 7320-37-8 | | HARMFUL / IRRITANT / CARCINOGENIC POTENTIAL | | * | 95,05,30,003 |
| EPOXYBISPHENOL RESIN | | | | | * | 89,09,04,011 |
| ERGOPOL 92 | | | | | * | 01,05,20,009/FDS |
| ERIOCHROME T BLACK | 1787-61-7 | | IRRITANT | | * | 95,06,21,002 |
| ERPAMINE 6 | 108-91-8 | | SEE CYCLOHEXYLAMINE | | | |
| ETHANOL | 64-17-5 | | FLAMMABLE | SOLVENT | * | |
| ETHANOL SODIUM SULFONATE 50% | | | SEE SODIUM HYDRO-2 ETHYL SULFONATE 50% | | | |
| ETHANOLAMINE | 141-43-5 | | SEE AMINO-2 ETHANOL | | | |
| ETHANOX | | | HARMFUL, IRRITANT | PHENOLIC ANTIOXIDANT | * | 090204007/MSDS |
| ETHER | 60-29-7 | | SEE DIETHYL ETHER | | | |
| ETHIDIUM BROMIDE | 1239-45-8 | | TOXIC | | * | 01,05,09,001 |
| ETHOMEEN TO 12 | | | CORROSIVE / HARMFUL | | * | |
| ETHOXY (20) CETYL ALCOHOL | 9004-95-9 | | HARMFUL | | * | 03,05,16,012/FDS |
| ETHYL 3-ETHOXYPROPIONATE | 763-69-9 | | IRRITANT | | * | 02,11,25,003/FDSWEB |
| ETHYL ACETATE | 141-78-6 | | FLAMMABLE | SOLVENT | * | |
| ETHYL ACRYLATE | 140-88-5 | | HARMFUL/FLAMMABLE | | * | |
| ETHYL ALCOHOL | 64-17-5 | | SEE ETHANOL | | | |
| ETHYL BROMOACETATE | 105-36-2 | | VERY TOXIC | | * | |
| ETHYL CHLOROFORMATE | 541-41-3 | | VERY TOXIC/CORROSIVE/REACTS WITH WATER | | * | 06,03,01,002/MSDS |
| ETHYLBENZEN | 100-41-4 | | VERY FLAMMABLE, HARMFUL | SOLVENT | * | 070919002-MSDS |
| ETHYLENE DIAMINE TETRAACETIC ACID | 60-00-4 | | IRRITANT | | * | 95,06,19,011 |
| ETHYLENE DIAMINE TETRAACETONITRILE | 5766-67-6 | | | | * | 01,01,10,009 |
| ETHYLENE GLYCOL DIMETHYL ETHER | | | SEE 1,2-DIMETHOXYETHANE | | | |
| ETHYLENE GLYCOL MONOETHYL ETHER | | | SEE 2-ETHOXYETHANOL | | | |
| ETHYLENE OXIDE | 75-21-8 | | TOXIC/IRRITANT/FLAMMABLE | | * | 95,06,01,022 |
| ETHYLENECHLORHYDRINE | 107-07-3 | | SEE 2-CHLOROETHANOL | | | |
| ETHYLENEDIAMINE | 107-15-3 | | CORROSIVE | BASE | * | 01,01,10,007 |
| ETHYLENEGLYCOL | 107-21-1 | | HARMFUL | | * | |
| ETHYLHEXYLAMINE | 104-75-6 | | TOXIC/CORROSIVE | AMINE | * | 070914002/MSDS |
| ETHYLMERCAPTAN | 75-08-1 | | HARMFUL, dangerous to the environment | SOLVENT | * | 090204006/MSDS |
| EUDRAGIT S12,5P | | | | | * | |
| FALSANE | 67485-29-4 | | HARMFUL | | * | |
| FEROCLEAN N-SF | | | | | * | 99,11,02,009 |
| FERROUS IRON | 7720-78-7 | | CORROSIVE | | * | 03,11,24,005/FDS |
| FILOXID | | | | INK | * | |
| FINOX DH N (PHOSPHORIC ACID / SOAP) | | | CORROSIVE | | * | 00,07,07,001 |
| FIRST PASS 5000 | | | IRRITANT | ACID | * | 100518002/MSDS |
| FLINT GROUP 1111 FLEXIPRINT COLOR RANGE | | | HARMFUL, IRRITANT | INK | * | 090317001/MSDS |
| FLUORFOLPET | | | IRRITANT | | * | 96,05,23,001 |
| FLUORINE | 7782-41-4 | | TOXIC/CORROSIVE | OXIDIZER | * | |
| FLUOROBORIC ACID | 16872-11-0 | | CORROSIVE/TOXIC | | * | 02,02,12,002 |
| FLYTANDE GRÖN | | | CORROSIVE | BASE | * | 00,02,21,001 |
| FORAPERLE | | | IRRITANT/HARMFUL/TOXIC | IRRITANT / NOXIOUS | * | 00,05,29,001 |
| FORMALDEHYDE | 50-00-0 | | SEE FORMOL | | | |
| FORMAMIDINE SULFINIC ACID | 1758-73-2 | | IRRITANT | | * | 080225001/MSDS |
| FORMIC ACID | 64-18-6 | | CORROSIVE | ACID | * | |
| FORMIC ALDEHYDE | 50-00-0 | | SEE FORMOL | | | |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | |
|--|------------|---|---------------------------|---|--------------------|
| FORMOL | 50-00-0 | TOXIC | SOLVENT | * | |
| FREONS | | NOT VERY TOXIC | CFC | * | |
| FUEL NAPHTHALINE | | | | * | |
| FUEL SOUTES | 68476-33-5 | TOXIC/FLAMMABLE | HYDROCARBURES MIX | * | 99,02,22,003 |
| FUELSOLV-FMG2900 | | | | * | 06.04.11.001/MSDS |
| FURFURAL | 98-01-1 | TOXIC/IRRITANT/CARCINOGEN | SOLVENT | * | 070912001/MSDS |
| FURFURYLAMINE | 617-89-0 | | BASE | * | 04,06,24,003/FDS |
| FURFURYLIC ALCOHOL | 98-00-0 | HARMFUL | SOLVENT | * | |
| FYRQUEL EHC (Rinçage) | | | | * | 04,06,23,001 |
| FYRQUEL EHC (Rinçage) | | | | * | 04,06,23,001 |
| GARDOBOND X 4650 | | IRRITANT | SURFACE PREPARATING AGENT | * | 081014001/MSDS |
| GAS OIL MARKER CONCENTRATE | | IRRITANT | SOLVENT | * | 090204004/MSDS |
| GASSER SOLUTION | | TOXIC,IRRITANT | BASE | * | 090325003/MSDS |
| GENAMIN LAP 100/D160Kg St-Drum,b (r.) | 5538-95-4 | CORROSIVE | BASE | * | 02,03,13,001 |
| GENAMINE TA 100 | 61790-33-8 | CORROSIVE/DANGEROUS TO THE ENVIRONMENT | FATTY AMINE | * | 98,10,27,003 |
| GEOLFLEX KOMP.B | | SEE DIPHENYLMETHANEDIISOCYANATE PREPOLYMERISED | PREPOLYMER | * | 100315003/MSDS |
| GERMANIUM TETRACHLORIDE | 10038-98-9 | IRRITANT | | * | 02,04,26,002 |
| GLUCONIC ACID | 526-95-4 | | WEAK ACID | * | |
| GLUTARAL | 111-30-8 | SEE GLUTARALDEHYDE | | | |
| GLUTARALDEHYDE | 111-30-8 | HARMFUL TO TOXIC / CORROSIVE/DANGEROUS TO THE ENVIRONMENT | | * | 95,06,20,002 |
| GLYCEROL EPICHLORHYDRINE | | SEE EPICHLORHYDRINE | | | |
| GLYCEROLGLYCID | 556-52-5 | SEE GLYCIDOL | | | |
| GLYCIDOL | 556-52-5 | CMR/TOXIC / HARMFUL / IRRITANT | | * | 95,05,29,018 |
| GLYCIDYLTRIMETHYLAMMONIUM CHLORIDE | 3033-77-0 | TOXIC/IRRITANT | | * | |
| GLYCINE | 56-40-6 | | WEAK ACID | * | |
| GLYCOLIC ACID | 79-14-1 | CORROSIVE/HARMFUL | CORROSIVE | * | 95,06,19,003 |
| GROTAN A | | TOXIC | | * | |
| HANUS SOLUTION | 7789-33-5 | FLAMMABLE / CORROSIVE | | * | 95,06,21,006 |
| HCFC 22 | 75-45-6 | | HCFC | * | 99,06,14,004 |
| HCR 840LIQUID (500) | | VERY TOXIC/CMR/CORROSIVE | | * | 061023001/MSDS |
| HELIANTHISE | 547-58-0 | SEE METHYL ORANGE | | | |
| HEPTYL ALCOHOL | | SEE HEPTANOL | | | |
| HERLI RAPID TW | | IRRITANT | | * | |
| HERLISIL | | CORROSIVE | | * | |
| HEXAFLUOROPHOSPHORIC ACID | 16940-81-1 | CORROSIVE/TOXIC | | * | 00,02,01,001 |
| HEXAMETHYLDISILOXANE | 107-46-0 | | | * | 06.07.04.002 /MSDS |
| HEXAMETHYLENE DIAMINE | 124-09-4 | CORROSIVE | BASE | * | 00,09,04,001 |
| HEXAMETHYLENETETRAMINE | 100-97-0 | HARMFUL/FLAMMABLE | BASE | * | |
| HEXANE | 110-54-3 | HARMFUL / FLAMMABLE | | * | 95,06,19,008 |
| HEXANEDIOIC ACID, POLYMER WITH N-(2-AMINOETHYL)-1,2-ETHANEDIAMINE AND (CHL | 25212-19-5 | IRRITANT | RESIN | * | 070208001/NEP |
| HEXYL ALCOHOL | 111-27-3 | SEE HEXANOL | | | |
| HFC 134a | 811-97-2 | | | * | 99,06,14,005 |
| HIPR | | IRRITANT | | * | 04,04,06,018/MSDS |
| HYDRATE SLURRY | | CORROSIVE | BASE | * | 01,10,15,001 |
| HYDRAZINE | 7803-57-8 | TOXIC | REDUCING AGENT/BASE | * | |
| HYDRAZINE HYDRATE 24% | 10217-52-4 | TOXIC/CORROSIVE/DANGEROUS TO THE ENVIRONMENT | CORROSIVE / NOXIOUS | * | 95,11,29,004 |
| HYDREX 200-800 | | | | * | 06.06.22.003/MSDS |
| HYDREX 4101 | | CORROSIVE | ACID/ POTASSIUM SALT | * | 070828001/MSDS |
| HYDRIODIC ACID | 10034-85-2 | | WEAK ACID | * | |
| HYDROBROMIC ACID | 10035-10-6 | CORROSIVE | | * | |
| HYDROCHLORIC ACID | 7647-01-0 | CORROSIVE | STRONG ACID | * | |
| HYDROFENE DISECO | | | | * | |
| HYDROFLUORIC ACID | 7664-39-3 | CORROSIVE / TOXIC | | * | 90,10,10,001 |
| HYDROFLUORIC ACID / AMMONIUM FLUORIDE | | CORROSIVE / TOXIC | ACID | * | 00,05,25,002 |
| HYDROFLUORIC ACID / NITRIC ACID | | CORROSIVE / TOXIC | ACID | * | 00,05,25,001 |
| HYDROFLUOSILIC ACID | 16961-83-4 | CORROSIVE | ACID | * | 02,11,29,001/FDS |
| HYDROGEN CYANAMIDE | 420-04-2 | TOXIC | | * | 01,11,06,001 |
| HYDROGEN CYANIDE | 74-90-8 | VERY TOXIC | ACID | * | 01,06,20,001 |
| HYDROGEN PEROXIDE | 7722-84-1 | CORROSIVE | OXIDIZER | * | |
| HYDROQUINONE | 123-31-9 | TOXIC | | * | |
| HYDROXYETHYLCELLULOSE | | | | * | 95,05,30,020 |
| HYPOPHOSPHOROUS ACID | 6303-21-5 | CORROSIVE | ACID, REDUCING AGENT | * | 091124001/MSDS |
| INDICATOR BUFFER TABLETS FOR THE DETERMINATION OF WATER HARDNESS WITH TITR | | IRRITANT/SENSITIZER | | * | 04,07,01,001 |
| INHIBITOR OP8402 | | IRRITANT | | * | 061009001/MSDS |
| INK 7190 | | IRRITANT | SOLVENT | * | |
| IODINE MONOCHLORIDE | 7790-99-0 | HARMFUL, CORROSIVE | ACID RELEASE | * | 080219001/MSDS |
| IODOMETHANE | 74-88-4 | TOXIC, IRRITANT, SENSITIZER | ALKYLATING AGENT | * | 100303001/MSDS |
| IRON | 7439-89-6 | FLAMMABLE | | * | 95,06,26,014 |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | | |
|---|------------|--|--------------------|---|--|-------------------|
| IRON (III) CHLORIDE ANHYDROUS | 7705-08-0 | CORROSIVE/HARMFUL/FLAMMABLE | | * | | |
| IRON (III) CHLORIDE HEXAHYDRATE | 10025-77-1 | CORROSIVE | | * | | |
| IRON PERCHLORATE | 13537-24-1 | TOXIC/FLAMMABLE | | * | | |
| IRON PERCHLORIDE | | SEE IRON TRICHLORIDE | | | | |
| IRON SULFATE HEPTAHYDRATE | 7782-63-0 | HARMFUL/IRRITANT | | * | | |
| IRON TRICHLORIDE ANHYDROUS | 7705-08-0 | CORROSIVE | | * | | 95,11,29,002 |
| IRON TRICHLORIDE HEXAHYDRATE | 10025-77-1 | CORROSIVE | | * | | 95,11,29,002 |
| IRON(II) CHLORIDE AND IRON(III) CHLORIDE IN HYDROCHLORIC ACID | | TOXIC/CORROSIVE/IRRITANT | | * | | |
| ISO FREE 977 LEAD FREE BASE | | CMR, IRRITANT | PAINT | * | | 080506001/MSDS |
| ISOAMYL ALCOHOL | 123-51-3 | SEE METHYL-3 BUTANOL | | | | |
| ISOBUTYL ACETATE | 110-19-0 | FLAMMABLE | SOLVENT | * | | |
| ISOBUTYLALCOHOL | 78-83-1 | SEE METHYL-2 PROPANOL | | | | |
| ISOBUTYLALDEHYDE | 78-84-2 | SEE 2-METHYLPROPANAL | | * | | |
| ISOCYANATES | | | | * | | |
| ISOPHORONE DIAMINE | 2855-13-2 | HARMFUL/CORROSIVE | | * | | |
| ISOPHORONE DIISOCYANATE | 4098-71-9 | TOXIC | | * | | 04,08,31,004/MSDS |
| ISOPROPANOL | 67-63-0 | SEE 2-PROPANOL | | | | |
| ISOPROPYL ACETATE | 108-21-4 | IRRITANT/FLAMMABLE | ESTER | * | | 98,08,17,005 |
| ISOPROPYL ALCOHOL | 67-63-0 | SEE 2-PROPANOL | | | | |
| ISOPROPYL TITANATE | 546-68-9 | IRRITANT/FLAMMABLE | | * | | 02,08,30,001 |
| JELLY FISH | | IRRITANT TO CORROSIVE | | * | | 05,09,27,001 |
| JET A-1 | | IRRITANT | | * | | 02,10,23,001/FDS |
| JKL (FeCl3 / Fe(SO4)3 / HCl) | | CORROSIVE | | * | | 02,02,19,001 |
| KATHON 893 MICROBIOCIDE | | TOXIC/CORROSIVE | | * | | |
| KATHON 930 BIOCID | | TOXIC/CORROSIVE | | * | | |
| KATHON WT14% BIOCID | | TOXIC/CORROSIVE | | * | | |
| KEBO ACIDE 710 | | CORROSIVE | | * | | 01,05,20,007/FDS |
| KEBOCOR N | | TOXIC BY INGESTION | | * | | 01,05,20,004/FDS |
| KEBOPLEX 102 | | | | * | | 01,05,20,001/FDS |
| KEBOPLEX III T | | IRRITANT | | * | | 01,05,20,005/FDS |
| KEBOSOL A | | | | * | | 01,05,20,006/FDS |
| KEBOSOL AN | | | | * | | 01,05,20,008/FDS |
| KEMAZUR 1215 | | | | * | | |
| KEMAZUR 1704 | | CORROSIVE | | * | | |
| KENWATER JKL | | CORROSIVE | | * | | 02,02,19,004 |
| KENWATER PAX-18 | | IRRITANT | | * | | 02,02,20,001 |
| KENWATER PAX-21 | 1327-41-9 | CORROSIVE | | * | | 02,02,19,003 |
| KILL NET FLOW | | | | * | | 01,05,20,010/FDS |
| KORASIT CC | | TOXIC/CMR/CORROSIVE/dangerous to the environment | ACID | * | | 070607001/MSDS |
| KRONOFLOC | | CORROSIVE | | * | | 02,11,26,001/MSDS |
| KYMENE | | W/OIC ACID, POLYMER WITH N-(2-AMINOETHYL)-1,2-ETHANEDIAMINE AND (CHLOROMET | | * | | |
| LACTIC ACID | 598-82-3 | IRRITANT | WEAK ACID | * | | |
| LANOLIC ACID | | FLAMMABLE | ACID | * | | 95,06,01,010 |
| LAURIC ACID | 143-07-7 | | WEAK ACID | * | | |
| LAURYL ALCOHOL | 112-53-8 | SEE DODECANOL | | | | |
| LAURYL BROMIDE | 143-15-7 | SEE 1-BROMODODECANE | | | | |
| LAVAZID 6027 | | IRRITANT | | * | | 01,11,06,002 |
| L-CYSTINE | 56-89-3 | | AMINED ACID | * | | 98,08,17,037 |
| LGA 9320 | | | VARNISH | * | | |
| LINCOMYCIN HYDROCHLORIDE | 859-18-7 | IRRITANT- POSSIBLE SENSITIZER | | * | | 04,12,07,001/FDS |
| LIQUID PARAFFIN | | IRRITANT | | * | | 04,12,07,008/MSDS |
| LITHIUM CARBONATE | 554-13-2 | HARMFUL | SALT | * | | |
| LITHIUM CHLORIDE | 7447-41-8 | IRRITANT, HARMFUL | SOLID | * | | 071004001-MSDS |
| LITHIUM HYDROXIDE | 1310-65-2 | CORROSIVE, Dangerous to the environment | BASE | * | | 071004002-MSDS |
| LITHIUM IN PARAFFIN OIL | 7429-93-2 | CORROSIVE | REDUCING AGENT | * | | 03,09,03,001/MSDS |
| LITHIUM TETRAHYDRIDOALUMINATE | 16853-85-3 | CORROSIVE,INFLAMMABLE,reacts with water | REDUCING AGENT | * | | 100927002/MSDS |
| L-LYSINE MONOHYDROCHLORIDE | 657-27-2 | | AMINED ACID | * | | 98,08,17,038 |
| LOCTITE 401 | | | CYANOACRYLATE GLUE | * | | 070625001/LET |
| LOCTITE 420 | | | GLUE | * | | 95,04,24,001 |
| LONZABAC | | HARMFUL, CORROSIVE | BIOCIDE, BASE | * | | 090910001/MSDS |
| L-SELECTRIDE SOLUTION | | SEE TRI-SEC-BUTYLHYDROBORATE | | * | | |
| LUCIDOL | | TOXIC/IRRITANT | | * | | |
| L'ZC 1051 | | | DILUTER | * | | |
| M GLASS ETCH | 1333-83-1 | SEE SODIUM BIFLUORIDE | ACID | * | | 100923001-b/FDS |
| M GLASS ETCH 3%/SULFURIC ACID 16% | | | ACID | * | | 100923001-a/FDS |
| MAGNESIUM CHLORIDE | 7786-30-3 | | | * | | 04,04,13,002/FDS |
| MAGNESIUM FLUOSILICATE | 16949-65-8 | TOXIC | | * | | 04,02,05,002/MSDS |
| MAGNESIUM SULFATE ANHYDROUS | 7758-98-7 | | | * | | 95,06,19,017 |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | | |
|--|------------|--|--|---|-------------------|-------------------|
| MAGNESIUM SULFATE HEPTAHYDRATE | 10034-99-8 | | | * | | 95,06,19,017 |
| MAGNESIUM SULFATE HYDRATE | 22189-08-8 | | | * | | 95,06,19,017 |
| MAGNUS 6-64 D | | | | * | | |
| MAGNUS 7-84 DISECO | | | | * | | |
| MAGNUS 92X3 | | | CORROSIVE | | BASE | 99,11,02,003 |
| MAKE UP SOLUTION | | | | * | BASE | 100722001/MSDS |
| MALEIC ACID | 110-16-7 | | HARMFUL | | ACID | |
| MALEIC ANHYDRIDE | 108-31-6 | | CORROSIVE/HARMFUL/IRRITANT/PROBABLY SENSITIZER | * | | 95,05,29,012 |
| MALIC ACID | 617-48-1 | | HARMFUL/IRRITANT | * | ACID | |
| MALONIC ACID | 141-82-2 | | IRRITANT | * | WEAK ACID | 98,08,17,013 |
| MANUGEL | | | IRRITANT | * | HYDROALCOOLIC GEL | 070518001/NEP |
| MARLOTHERM SM | 26898-17-9 | | | * | | 02,12,10,001/FDS |
| MBTCI | 1118-46-3 | | SEE n-BUTYL TIN TRICHLORIDE | * | | |
| MCA 80% IN METHANOL | | | TOXIC | * | | 01,10,08,001 |
| MCA 80% IN WATER | | | TOXIC | * | | 01,10,08,002 |
| MDI | 101-68-8 | | SEE METHYLENE DI-p-PHENYL DIISOCYANATE | | | |
| MEGAPOSIT (R) SPR 3012 PHOTORESIST | | | | * | | 04,06,24,001/MSDS |
| MELDRUM'S ACID | 2033-24-1 | | IRRITANT | * | ACID | 98,08,17,011 |
| MERCURIC ACETATE | 1600-27-7 | | VERY TOXIC | * | | 95,08,04,001 |
| MERCURY | 7439-97-6 | | TOXIC | * | | |
| METHACRYLIC ACID | 79-41-4 | | CORROSIVE | * | ACID | 95,05,30,030 |
| METHACRYLONITRILE | 126-98-7 | | TOXIC/PROBABLY SENSITIZER/HIGHLY FLAMMABLE | * | | 95,06,21,007 |
| METHANESULFONIC ACID | 75-75-2 | | CORROSIVE | * | ACID | |
| METHANOL | 67-56-1 | | TOXIC/FLAMMABLE | * | SOLVENT | |
| METHYL METHACRYLATE | 80-62-6 | | IRRITANT/FLAMMABLE | * | | 95,05,30,021 |
| METHYL 10-UNDECENOATE | 1731-86-8 | | SEE METHYL UNDECYLENATE | | | |
| METHYL ACETATE | 79-20-9 | | FLAMMABLE | * | ESTER | 95,06,19,001 |
| METHYL BROMIDE | | | SEE BROMOMETHANE | | | |
| METHYL FORMATE | 107-31-3 | | HARMFUL, IRRITANT | * | | 071003002-MSDS |
| METHYL ISOBUTYL CETONE | 108-10-1 | | SEE 4-METHYL-2-PENTANONE | | | |
| METHYL ISOCYANATE | 624-83-9 | | IRRITANT | * | REACTS WITH WATER | |
| METHYL ORANGE | 547-58-0 | | TOXIC | * | | |
| METHYL PARABEN | 99-76-3 | | SEE PARAHYDROXYBENZOATE DE METHYLE | | | |
| METHYL PARAHYDROXYBENZOATE | 99-76-3 | | | * | | 95,05,29,003 |
| METHYL TERTIOBTYL ETHER | 1634-04-4 | | FLAMMABLE/IRRITANT | * | | 95,05,30,036 |
| METHYL UNDECYLENATE | 1731-86-8 | | WEAK IRRITANT | * | | 95,06,02,002 |
| METHYL VINYLL KETONE | 78-94-4 | | SEE BUTENONE | * | | |
| METHYLACRYLATE | 96-33-3 | | HARMFUL/IRRITANT | * | | 070927003/MSDS |
| METHYLAMINE, 40 wt.% solution in water | 74-89-5 | | HARMFUL/CORROSIVE/FLAMMABLE | * | | |
| METHYLCHLOROACETATE | 96-34-4 | | TOXIC | * | ESTER | 01,10,08,001 |
| METHYLENE BLUE | 61-73-4 | | HARMFUL | * | | 95,06,21,001 |
| METHYLENE CHLORIDE | 75-09-2 | | SEE DICHLOROMETHANE | | | |
| METHYLENE IODIDE | 75-11-6 | | HARMFUL, IRRITANT | * | SOLVENT | 081112001/MSDS |
| METHYLENE DI-p-PHENYL DIISOCYANATE | 101-68-8 | | HARMFUL | * | NOXIOUS | |
| METHYLETHYL KETONE | 78-93-3 | | SEE 2-BUTANONE | | | |
| METHYLETHYLKETONE PEROXIDE | 1338-23-4 | | TOXIC/FLAMMABLE | * | STRONG OXIDIZER | 95,03,09,001 |
| METHYL-ISOTHIAZOLINONE | 2682-20-4 | | CORROSIVE | * | | 01,10,04,001 |
| MEYER-NITRO-UNIVERSAL-VERDÜNNUNG ENGL | | | IRRITANT | * | | 03,02,25,002/FDS |
| MIBC | 108-11-2 | | SEE 4-METHYLPENTAN-2-OL | * | | |
| MICROPOSIT (R) 351 DEVELOPER | | | | * | | 04,06,24,001/MSDS |
| MICROPOSIT (R) PRIMER | 999-97-3 | | CORROSIVE | * | | 03,03,31,003/FDS |
| MICROPOSIT EC SOLVENT | | | IRRITANT | * | | 03,03,31,005/FDS |
| MICROPOSIT MF-321 DEVELOPER | | | | * | | 04,06,24,001/MSDS |
| MICROPOSIT S1818 PHOTO RESIST | 108-65-6 | | IRRITANT | * | | 03,03,31,002/FDS |
| MICROTREAT 2700 | | | HARMFUL, IRRITANT, Toxic to the environment | * | BIOCIDE | 080618001/MSDS |
| MIDAS 98.2 | | | SEE IODOMETHANE | * | | 100303001/MSDS |
| MILLISIL C6 | | | | * | | |
| MIXED ACID ETCH 1:1:2 | | | CORROSIVE | * | ACID | |
| MIXED ACID ETCH 10:1:4:1 | | | TOXIC/CORROSIVE | * | ACID | |
| MIXED ACID ETCH 3:1:1 | | | TOXIC/CORROSIVE | * | ACID | |
| MIXED ACID ETCHANT 4:2:2 | | | TOXIC/CORROSIVE | * | ACID | |
| MIXTURE OF METHYLISOTHIAZOLONE AND CHLOROMETHYLISOTHIAZOLONE | 55965-84-9 | | TOXIC, CORROSIVE, SENSITIZER | * | BIOCIDE | 080225002/MSDS |
| MOBIL DTE 24 | | | | * | HYDRAULIC FLUID | 080516001/MSDS |
| MODIFIED COPRAH FAT ACID | | | IRRITANT | * | | 95,06,01,001 |
| MOLECULAR SIEVES 10A | | | IRRITANT | * | | 04,04,06,010/MSDS |
| MONOBUTYL TIN TRICHLORIDE | 1118-46-3 | | CORROSIVE | * | | 02,08,30,003 |
| MONOISOPROPYLAMINE | 75-31-0 | | SEE 1-AMINO-2-PROPANOL | | | |
| MONOMETHYLHYDRAZINE ANHYDROUS | 60-34-4 | | TOXIC/CORROSIVE | * | REDUCING AGENT | 02,08,01,001 |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | |
|--|-------------|---|--------------------------------|---|-----------------------------|
| MORPHOLINE | 110-91-8 | CORROSIVE/HARMFUL/FLAMMABLE | | * | |
| MTBE | 1634-04-4 | SEE METHYL TERTIOBUTYL ETHER | | | |
| MYRISTIC ACID | 544-63-8 | IRRITANT | ACID | * | 95,05,29,005 |
| N-(2-AMINOETHYL) ETHANOLAMINE | 111-41-1 | SEE N-(2-HYDROXYETHYL) ETHYLENEDIAMINE | | | |
| N-(2-HYDROXYETHYL) ETHYLENEDIAMINE | 111-41-1 | CORROSIVE | BASE | * | 95,05,30,004 |
| N-(FLUORODICHLOROMETHYLTHIO)-PHTHALIMIDE | | IRRITANT | | * | |
| N,N-DIBUTYLFORMAMIDE | 761-65-9 | IRRITANT | SOLVENT | * | 071003001-MSDS |
| N,N'-DICYCLOHEXYLCARBODIIMIDE | 538-75-0 | SEE DCC | | | |
| N,N'-DIETHYL-1,3-DIAMINOPROPANE | 104-78-9 | SEE 3-DIETHYLAMINO-1-PROPYLAMINE | | | |
| N,N'-DIETHYLAMINO-2-ETHANOL | 100-37-8 | IRRITANT | | * | 95,05,29,007 |
| N,N'-DIETHYL-N,N'-DIMETHYLPROPANE DIAMINE | | HARMFUL / CORROSIVE | | * | 95,05,30,009 |
| N,N-DIMETHYLACETAMIDE | | | | * | 06.03.01.001/MSDS |
| N,N-DIMETHYLDODECYLAMINE | 112-18-5 | HARMFUL / CORROSIVE | | * | 95,06,01,005 |
| N,N'-ETHYLENETHIOUREA | 96-45-7 | TOXIC | BASE | * | 98,08,17,033 |
| NALCO 4360 TAB | | HARMFUL,CORROSIVE,SENSIBILISANT, dangerous to the environment | BIOCIDE | * | 090819001/FDS |
| NALCO 71101 | | SENSITIZER/DANGEROUS TO THE ENVIRONMENT | anti-foaming agent | * | 070813002/MSDS |
| NALCO 7208 | | CORROSIVE | BASE | * | 070810002/MSDS |
| NALCO 72100 | | IRRITANT/SENSITIZER | REDUCING AGENT | * | 070813001/MSDS |
| NALCO 73190 | | SEE Sodium 4(or 5)-methyl-1H-benzotriazolide | BASE & REDUCING AGENT | * | 080229001/MSDS |
| NALCO 74724 | | | | * | 06.03.29.001/FDS |
| NALCO 74836 | | | | * | 06.03.29.002/FDS |
| NALCO 77352 | | CORROSIVE, SENSIBILISANT | BIOCIDE,ACID | * | 090819002/FDS |
| NALCO 77393 | | | BIOLOGICAL DETERGENT | * | 090819003/FDS |
| NALCO DUSTFOAM PLUS | | IRRITANT, Dangerous to the environment | ACID | * | 080808001/MSDS |
| NALCO® 7320 | | CORROSIVE | ACID | * | |
| N-ALKYL(C12-C22)-1,3-DIAMINOPROPANE | 90640-44-1 | CORROSIVE/DANGEROUS TO THE ENVIRONMENT | OSIVE/DANGEROUS TO THE ENVIRON | * | 98,11,09,001 |
| NANO TM MMA (L7,5) MAA COPOLYMER SERIES RESISTS | | | | * | 04,06,24,001/MSDS |
| NAPHTHA(petroleum) | 64742-48-9 | HARMFUL | | * | 06,02,03,003/MSDS |
| NAPHTHYLAMINE | 134-32-7 | TOXIC/DANGEROUS TO THE ENVIRONMENT | | * | 95,06,21,008 |
| N-BROMOSUCCINIMIDE | 128-08-5 | HARMFUL/IRRITANT | | * | |
| n-BUTYLLITHIUM | 109-72-6 | | STRONG REDUCING AGENT | * | 071004003-MSDS |
| n-BUTYLLITHIUM in HEXANE | | SEE n-BUTYLLITHIUM | | * | 071004003-MSDS |
| n-BUTYLITIN TRICHLORIDE | 1118-46-3 | CORROSIVE,dangerous to the environment | ACID, ALKYLATING AGENT | * | 090430002/MSDS |
| N-COCOALKYL PROPYLENE DIAMINE | | CORROSIVE | BASE | * | 95,06,01,006 |
| NDM | | SEE NORMAL DODECYLMERCAPTAN | | | |
| NEOMYCIN SULFATE | 1405-10-3 | SENSITIZER | | * | 04,12,07,004/MSDS |
| NEOPENTYLGLYCOL | | SEE 2,2-DIMETHYL-1,3-PROPANEDIOL | | | |
| NEOSTANN U-220H | | SEE Dibutylbis(pentane-2,4-dionato-O,O')tin | | * | |
| NET INOX N | | SEE HYDROFLUORIC ACID / NITRIC ACID | | | |
| NETTOR AL 12 (HF/H2SO4) | | CORROSIVE/TOXIC | ACID | * | 01,09,28,001 |
| NEUTRACIT -2 CAUSTIC NEUTRALIZER | | IRRITANT | | * | 04,04,06,011/MSDS |
| n-HEPTANE | 142-82-5 | FLAMMABLE | | * | 95,05,30,033 |
| NICKEL CHLORIDE HEXAHYDRATE | 7791-20-0 | TOXIC | | * | 01,10,23,001 |
| NICKEL SALTS | | IRRITANT OR TOXIC | | * | |
| NINHYDRINE | 485-47-2 | HARMFUL/IRRITANT | | * | 02,08,29,005 |
| NITRIC ACID | 7697-37-2 | CORROSIVE | STRONG ACID | * | |
| NITRIC ACID 15 % AND HYDROFLUORIC ACID | | CORROSIVE / TOXIC | ACID | * | 96,01,24,001 |
| NITRIC ACID ETCHANT SF 68-01 VLSI SELECTIPUR | | | | * | 04,06,24,002/MSDS |
| NITROBENZENE | 98-95-3 | TOXIC | SOLVENT | * | 070423001/MSDS |
| NITROGENE DIOXYDE | 10102-44-00 | IRRITANT | | * | 00,06,06,001 |
| NITROMETHANE | 75-52-5 | HARMFUL | | * | 98,08,17,044 |
| N-METHYLPYRROLIDONE | 872-50-4 | IRRITANT | | * | 98,08,17,045 |
| n-OCTANAL | 124-13-0 | | | * | 03,05,16,011/FDS |
| NONANAL | 124-19-6 | IRRITANT | | * | 03,05,16,006/FDS |
| NONIVAMIDE | | IRRITANT (PAVA SPRAY or GAS) | | * | 05,06,23,001/CAS |
| NORMAL DODECYLMERCAPTAN | 112-55-0 | IRRITANT | | * | 05,07,08,001/FDS |
| NORPOL FI-165 | | | | * | 06.06.22.002/MSDS |
| NORPOL PEROXIDE 18 | | | | * | 06.06.22.001/MSDS |
| N-VINYL-2-PYRROLIDINONE | 88-12-0 | HARMFUL | | * | 99,02,03,002 |
| O-(BENZOTRIAZOL-1-YL)-N,N,N',N'-TETRAMETHYLURONIUM TETRAFLUOROBORATE | | SEE TBTU | | | |
| OC SPRAY or GAS | | IRRITANT | GAS | * | 05,08,30,001-2/biblio-tests |
| OCTADECANOIC ACID | 57-11-4 | SEE STEARIC ACID | | | |
| OCTADECYL ALCOHOL | 112-92-5 | SEE 1-OCTADECANOL | | | |
| OCTANOYL CHLORIDE | 111-64-8 | CORROSIVE | CHLORURE D ACIDE | * | 98,08,17,023 |
| OCTYL ALCOHOL | 111-87-5 | SEE OCTANOL | | | |
| OCTYLPHOSPHONIC ACID | 4724-48-5 | FLAMMABLE, TOXIC, CORROSIVE | ACID | * | 091211001/NFP |
| OLEIC ACID | 112-80-1 | IRRITANT | | * | 95,05,29,002 |
| OLEIC AND TALLOWETHOXYL AMINE | | CORROSIVE / HARMFUL | | * | 95,06,02,001 |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | |
|--|--------------|---|----------------------------------|---|-----------------------------|
| OLEORESIN CAPSICUM | | SEE OC SPRAY | | | |
| OLEOYL ALCOHOL | 143-28-2 | IRRITANT | | * | 95,05,30,013 |
| OLEOYL CHLORIDE | 112-77-6 | IRRITANT | | * | 95,05,30,037 |
| OLEUM 20% | 98,10,27,001 | CORROSIVE | STRONG ACID | * | 98,10,27,001 |
| OPALKLEEN DISECO | | | | * | |
| OPD 5262 | 75-59-2 | HARMFUL/IRRITANT | | * | 03,03,31,004/FDS |
| O-TOLUENEDIAMINE-LAA | | | MANUFACTURING OF URETHANNES ACID | * | 05,07,06,003/MSDS |
| OXALIC ACID | 144-62-7 | HARMFUL | | * | |
| OXIRANE-2-METHANOL | 556-52-5 | SEE GLYCIDOL | | | |
| OXYTETRACYCLINE DIHYDRATE | | HARMFUL/IRRITANT | | * | 04,12,07,005/MSDS |
| P3 FERROLIN | | CORROSIVE | BASE | * | 070718001/FDS |
| P3 PREVOX 675 | | | | * | 99,11,02,008 |
| PAINTEX 1967 | | CORROSIVE, HARMFUL | | * | 071120001-MSDS |
| PALLA ASSIST | | SEE PALLADIUM (II) TETRAAMINE CHLORIDE | | * | |
| PALLADIUM (II) TETRAAMINE CHLORIDE | 13815-17-3 | IRRITANT, SENSITIZER | BASE | * | 100818001/MSDS |
| PALMITOYL CHLORIDE | 112-67-4 | CORROSIVE | CHLORURE D ACIDE | * | 98,08,17,020 |
| para-ANISYL ALCOHOL | 105-13-5 | SEE P-METHOXYBENZOIC ALCOHOL | | | |
| PARAHYDROXYBENZALDEHYDE | 123-08-0 | IRRITANT | | * | 98,08,17,048 |
| PARA-TOLUENE SULFONIC ACID | 6192-52-5 | IRRITANT | | * | 95,05,30,006 |
| PARCO 2074 | | CORROSIVE | BASE | * | 99,11,02,001 |
| PARCOLUBRITE 81 | | | | * | 99,11,02,002 |
| PÂTE DECAPANTE PX100 | | TOXIC/CORROSIVE | ACID | * | |
| PAVA SPRAY or GAS | | IRRITANT | | * | 05,08,30,001-2/biblio-tests |
| PAX | 2720-73-2 | SEE POTASSIUM AMYL XANTHATE | | * | |
| PDS PH-1025 PHOSPHORUS PLANAR DIFFUSION SOURCE | | IRRITANT | | * | 04,04,06,012/MSDS |
| PENNGAR MANUCHLORE 500 | | CORROSIVE, dangerous to the environnement | BASE OXYDIZER | * | 101123001/FDS |
| PENTANEDIAL | 111-30-8 | SEE GLYTARALDEHYDE | | | |
| PENTENAL | 107-02-8 | SEE ACRYLALDEHYDE | | * | |
| PENTYL ACETATE | 628-63-7 | SEE AMYL ACETATE | | | |
| PEPITE NEUTRE | | FLAMMABLE, IRRITANT, SENSITIZER | ACRYLATES | * | 080326002/FDS |
| PEPPER SPRAY or GAS | | IRRITANT | IRRITANT | * | 05,08,30,001-2/biblio-tests |
| PERACETIC ACID | 79-21-0 | CORROSIVE/HARMFUL/FLAMMABLE | | * | |
| PERCHLOROETHYLENE | 127-18-4 | SEE TETRACHLOROETHYLENE | | | |
| PERCHLORIC ACID | 7601-90-3 | CORROSIVE/FLAMMABLE | ACID | * | |
| PERCY 141 | | CORROSIVE | BASE | * | 99,06,14,002 |
| PERFLUOROPROPANE | 76-19-7 | | | * | 04,04,06,008/MSDS |
| PERGAQUICK C 24 AP | | TOXIC | | * | 01,01,23,009 |
| PERIODIC ACID | 1045-60-9 | CORROSIVE/OXIDIZING | WEAK ACID | * | |
| PEROXANE ME 50L | | CORROSIVE (SEE METHYLETHYLKETONPEROXIDE AND DIMETHYLPHALAT) | | * | 01,01,23,008 |
| PEROXYBOOST P420 | | | | * | 02,03,04,001 |
| PESTICIDES | | | | * | 03,09,16,003/DOSSIER |
| PETROLEUM ETHER | | HARMFUL/FLAMMABLE | | * | 95,06,21,005 |
| PHENETYL ALCOHOL | 60-12-8 | HARMFUL/IRRITANT | ALCOHOL | * | 95,06,21,003 |
| PHENOL RED | 143-74-8 | IRRITANT | | * | 91,01,03,001 |
| PHENOLPHTHALEINE | 77-09-8 | FLAMMABLE | | * | 03,08,29,001/FDS |
| PHENYLTRICHLOROSILANE | 98-13-5 | TOXIC/CORROSIVE | | * | 04,08,31,003/MSDS |
| PHOSGENE | 75-44-5 | TOXIC/IRRITANT | | * | 070709001/NEP |
| PHOSPHINE | 7803-51-2 | FLAMMABLE - TOXIC BY INHALATION | | * | 00,05,25,002 |
| PHOSPHORIC ACID | 7664-38-2 | CORROSIVE | STRONG ACID | * | 98,08,17,046 |
| PHOSPHORIC ACID / AMMONIUM FLUORIDE | | CORROSIVE/TOXIC | CORROSIVE/TOXIC | * | 03,05,12,001/FDS |
| PHOSPHORUS OXIDE CHLORIDE | 10025-87-3 | CORROSIVE | | * | |
| PHOSPHORUS PENTAOXIDE | 1314-56-3 | CORROSIVE | | * | 03,09,16,002/MSDS |
| PHOSPHORUS TRICHLORIDE | 2125683 | CORROSIVE | | * | 98,08,17,046 |
| PHOSPHORUS TRIFLUORIDE | 7783-55-3 | TOXIC/CORROSIVE | | * | 99,11,15,001 |
| PHOSPHORYL CHLORIDE | 10025-87-3 | VERY TOXIC/ORROSIVE/REACTS WITH WATER | | * | 95,05,29,020 |
| PHTHALIC ANHYDROUS | 85-44-9 | IRRITANT | | * | 98,08,17,049 |
| PIPERAZINE ANHYDROUS | 110-85-0 | CORROSIVE | | * | 03,03,31,001/FDS |
| PIPERIDINE | 110-89-4 | TOXIC/FLAMMABLE | | * | 100315003/MSDS |
| PIX-1400 | 872-50-4 | IRRITANT | | * | |
| pMDI | | SEE DIPHENYLMETHANE DIISOCYANATE PREPOLYMERISED | PREPOLYMER | * | |
| P-METHOXYBENZYL ALCOHOL | 105-13-5 | HARMFUL | ALCOHOL | * | 99,11,02,006 |
| POLICLEAN 251 | | CORROSIVE | | * | |
| POLY HEXAMETHYLENE DIISOCYANATE | | SEE DESMODUR N-3200 | | | |
| POLYALUMINIUM CHLORIDE | 1327-41-9 | CORROSIVE | ACID SOLUTION | * | 070914001/MSDS |
| POLYALUMINIUMHYDROXIDECHLORIDESULFATE | 39290-78-3 | IRRITANT | | * | 05,07,06,001/MSDS |
| polyamidoamine-epichlorhydrine (PAAE) resin in water | | IRRITANT | RESIN | * | 080225003/MSDS |
| POLYCAT 5 | | SEE BIS(2-DIMETHYLAMINOETHYL)(METHYL)AMINE | | * | |
| POLYETHYLENEGLYCOLS | 25322-68-3 | | | * | |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | | |
|---|------------|---|-------------------|---|---|----------------------------|
| POLYOLS | | | | * | | |
| POLYURETHANNE | | | | * | | |
| POLYVINYL ALCOHOL | 9002-89-5 | | | * | | |
| PORTLAND CEMENT | 65997-15-1 | FLAMMABLE - IRRITANT | | * | | 04,08,31,001/MSDS |
| POTASH | 1310-58-3 | CORROSIVE | STRONG BASE | * | | |
| POTASH DETERGENT | | CORROSIVE | | * | | |
| POTASSIUM AMYL XANTHATE | 2720-73-2 | HARMFUL,IRRITANT,REACTS WITH WATER | BASE | * | | 090616001/MSDS |
| POTASSIUM BICARBONATE | 298-14-6 | | POTASSIUM SALT | * | | |
| POTASSIUM BICHROMATE | 7778-50-9 | VERY TOXIC/CORROSIVE/CMR/DANGEROUS TO THE ENVIRONMENT | POTASSIUM SALT | * | | |
| POTASSIUM BIFLUORIDE | 7789-29-9 | TOXIC/CORROSIVE | FLUORIDE | * | * | 04,02,05,001/MSDS-Prometra |
| POTASSIUM BISULFITE | 16731-55-8 | IRRITANT | POTASSIUM SALT | * | | |
| POTASSIUM CARBONATE ANHYDROUS | 584-08-7 | HARMFUL/IRRITANT | POTASSIUM SALT | * | | |
| POTASSIUM CHLORIDE | 7789-00-6 | IRRITANT | POTASSIUM SALT | * | | |
| POTASSIUM CHROMATE | 7789-00-6 | TOXIC/IRRITANT/DANGEROUS TO THE ENVIRONMENT | OXIDIZER | * | | |
| POTASSIUM DICHROMATE | 7778-50-9 | SEE POTASSIUM BICHROMATE | | | | |
| POTASSIUM DIHYDROGEN PHOSPHATE | 7778-77-0 | | | * | | |
| POTASSIUM DIPHOSPHATE | 7320-34-5 | | SOLIDE | * | | 06,02,02,002/MSDS |
| POTASSIUM DISULFITE | 16731-55-8 | IRRITANT | POTASSIUM SALT | * | | |
| POTASSIUM FLUORIDE | 7789-23-3 | TOXIC | | * | | 04,08,31,002/FDS |
| POTASSIUM HEXACYANOFERRATE(II) TRIHYDRATE | 14459-95-1 | IRRITANT | | * | | 02,02,12,001 |
| POTASSIUM HEXACYANOFERRATE(III) | 13746-66-2 | IRRITANT | | * | | 02,02,12,001 |
| POTASSIUM HEXAFLUORIDE ANTIMONATE | 16893-92-8 | CORROSIVE/IRRITANT | | * | | |
| POTASSIUM HEXAFLUORIDE PHOSPHATE | 17084-13-8 | IRRITANT | | * | | |
| POTASSIUM HEXAFLUORIDE TITANATE | 16919-27-0 | IRRITANT | | * | | |
| POTASSIUM HYDROGEN PHTHALATE | 877-24-7 | | | * | | |
| POTASSIUM HYDROXIDE | 1310-58-3 | SEE POTASH | | | | |
| POTASSIUM IODATE | 7758-05-6 | HARMFUL/IRRITANT/OXIDIZING | | * | | 95,06,26,004 |
| POTASSIUM IODIDE | 7681-11-0 | | POTASSIUM SALT | * | | 98,08,17,036 |
| POTASSIUM METABISULFITE | 16731-55-8 | SEE POTASSIUM DISULFITE | | | | |
| POTASSIUM NITRATE | 7757-79-1 | OXIDIZING | POTASSIUM SALT | * | | |
| potassium O-pentyl dithiocarbonate | 2720-73-2 | SEE POTASSIUM AMYL XANTHATE | | * | | |
| POTASSIUM PERMANGANATE | 7722-64-7 | HARMFUL/OXIDIZING | OXIDIZER | * | | |
| POTASSIUM PERSULFATE | 7727-21-1 | HARMFUL/OXIDIZING | POTASSIUM SALT | * | | |
| POTASSIUM SILICATE | | | POTASSIUM SALT | * | | |
| POTASSIUM SULFITE | 10117-38-1 | IRRITANT | | * | | |
| POTASSIUM tert-BUTOXIDE | 865-47-4 | CORROSIVE/FLAMMABLE | VERY STRONG BASE | * | | 99,06,02,001 |
| PQ100 | | | | * | | 02,03,04,001 |
| PRAESTOL 2640 | | | ACRYLIC COPOLYMER | * | | 99,09,20,001 |
| PRAESTOL 852 BC | 75150-29-7 | MINIUM,N,N,N-TRIMETHYL-3-[(1-OXO-2-PROPENYL)AMINO],CHLORIDE, POLYMER WITH | | * | | |
| PREPARKER 43A | | | | * | | 99,11,02,007 |
| PREPARKER 43B | | | | * | | 99,11,02,004 |
| PREVENTOL CMK | 59-50-7 | TOXIC/IRRITANT | | * | | |
| PREVENTOL GDA/IT | | CORROSIVE,SENSITIZER | BIOCIDE | * | | 091116002/MSDS |
| PRIMASEPT MED | | IRRITANT /FLAMMABLE | | * | | 06,06,14,001/MSDS |
| PROCAP AV | | TOXIC/CORROSIVE | ACID | * | * | |
| PROCIV 300 | | CORROSIVE | ACID | * | | |
| PROLINK MOSQUITO GROWTH FACTORS | | HARMFUL | INSECTICIDE | * | | 100513001/MSDS |
| PROLINK PELLETS MOSQUITO GROWTH FACTORS | | HARMFUL | INSOLUBLE SOLID | * | | 100513002/MSDS |
| PROLINK PROSAND MOSQUITO GROWTH FACTORS | | HARMFUL | INSOLUBLE SOLID | * | | 100513003/MSDS |
| PROLINK XR BRIQUETS MOSQUITO GROWTH FACTORS | | HARMFUL | INSOLUBLE SOLID | * | | 100513004/MSDS |
| PROPANOIC ACID | 79-09-4 | CORROSIVE | CARBOXYLIC ACID | * | | 070919001/MSDS |
| PROPARGYL ALCOHOL | 107-19-7 | TOXIC/FLAMMABLE | ALCOHOL | * | | |
| PROPERASE L 1600 | | HARMFUL/IRRITANT | | * | | |
| PROPIONIC ACID | 79-09-4 | SEE PROPANOIC ACID | | * | | |
| PROPYL ALCOHOL | 71-23-8 | SEE PROPANOL-1 | | * | | |
| PROPYLENE OXIDE | 75-56-9 | CMR, TOXIC, IRRITANT | ALKYLATING AGENT | * | | 100114001/MSDS |
| PROPYLENEGLYCOL | 57-55-6 | SEE PROPANEDIOL | | * | | 98,08,07,041 |
| PROPYLENEGLYCOL MONOETHYLETER | 107-98-2 | SEE 1-METHOXY2-PROPANOL | | * | | |
| PROSOLV 500 | | IRRITANT | | * | | |
| PRÜFLÖSUNG 10 | 7789-23-3 | TOXIC | | * | * | 02,10,02,001/FDS |
| PSE 300 | | | | * | | 04,04,13,003/FDS |
| P-TOLUENESULFONYL ISOCYANATE | 4083-64-1 | IRRITANT/REACTS WITH WATER | | * | | 06,08,17,001/EXP |
| PUMICE STONE | | | | * | | 95,06,26,010 |
| PYRALENE | | | | * | | |
| PYRIDINE | 110-86-1 | HARMFUL/FLAMMABLE | BASE | * | | |
| QUICKLIME | 1305-78-8 | SEE CALCIUM OXIDE | | * | | |
| QUINOLINE | 91-22-5 | CMR, HARMFUL,IRRITANT | BASE | * | | 090204002/MSDS |
| R-507 | | | | * | | 99,06,14,007 |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | |
|---|------------|--|-----------------------|---|-------------------|
| RADIAMINE G161 | | CORROSIVE | | * | |
| RAFANOL P 80/85 | | | | * | |
| RECHARGEABLE LITHIUM BATTERY | | IRRITANT, SENSITIZER | | * | 080708002/MSDS |
| RECHARGEABLE POLYMER LI-ION BATTERY | | CMR, IRRITANT | | * | 080709001/MSDS |
| REDICOTE N39L | | CORROSIVE | BASE | * | 99,03,22,001 |
| REFINED COPRAH OIL | | | | * | 95,06,01,013 |
| RESINE 3052 JEANNEAU | | | | * | |
| RETICULATE THERMORIGID POLYESTER | | | | * | |
| RHODAMEEN OS 12 | | CORROSIVE / HARMFUL | | * | |
| RIDELBAX MILAN 1 RH 848 S1 | | IRRITANT | PAINT | * | 080326007/FDS |
| ROUND UP BIOACTIVE | | IRRITANT | ACID | * | 100512002/MSDS |
| SACHTOKLAR | | SEE POLYALUMINIUMHYDROXIDECHLORIDESULFATE | | | |
| SALICYLIC ACID | 69-72-7 | IRRITANT | WEAK ACID | * | 98,08,17,014 |
| SALPETERSAURE-Atzmischung SF 70,5-1,5 STAUBKLASSE 0-2 MOS | | | | * | 04,06,24,001/MSDS |
| SC-2000 CEMENT | | CMR, IRRITANT | SOLVENT | * | 100217002/MSDS |
| SCANNER | | HARMFUL | | * | 05,01,20,013/FDS |
| SEDIPUR CF 5001 | | | FLOCCULANT AGENT | * | 061222001/MSDS |
| SEFTOL E.F, 95 REYDEL | | | | * | |
| SELENIC ACID | 7783-00-8 | TOXIC, CORROSIVE | OXIDIZER | * | 100827001/MSDS |
| SELENIUM | 7782-49-2 | TOXIC | | * | 95,06,26,011 |
| SHELL MWS ADDITIVE B201 | | CORROSIVE | | * | 01,12,04,001 |
| SHELL TELLUS 32 | | IRRITANT | | * | |
| SHUREGEL 5 | | CMR, CORROSIVE | BASE | * | 100722002/MSDS |
| SIGMACOTE | | | SOLVENT | * | 04,04,06,014/MSDS |
| SILICAGEL | | | | * | |
| SILICIUM DIOXIDE | 14808-60-7 | | | * | 06,01,17,001/MSDS |
| SILICIUM TETRACHLORIDE | 10026-04-7 | IRRITANT | | * | 02,04,26,001/NEP |
| SILICON TETRAFLUORIDE | 07783-61-1 | CORROSIVE/TOXIC | | * | 02,07,11,001 |
| SILVER CONDUCTIVE PASTE | | HARMFUL, IRRITANT | | * | 100816002/MSDS |
| SILVER NITRATE DILUATED | | IRRITANT | | * | 95,08,04,002 |
| SILVER NITRATE PURE | 7761-88-8 | CORROSIVE | | * | 95,08,04,002 |
| SIMULSOL 58 | 9004-95-9 | HARMFUL | | * | 03,05,16,004/FDS |
| SINESTO B | | HARMFUL, CORROSIVE, Dangerous to the environment | BASE | * | 080821001/MSDS |
| SK10160 | 64742-48-9 | HARMFUL | | * | 06,02,003/MSDS |
| SKYDROL | | IRRITANT | | * | 95,12,13,001 |
| SLAKED LIME | 1305-62-0 | SEE CALCIUM HYDROXIDE | | * | |
| SODA | 1310-73-2 | CORROSIVE | STRONG BASE | * | 90,09,24,002 |
| SODA DETERGENT | | CORROSIVE | | * | |
| SODA METHYLATE 30 % IN MeOH | | TOXIC / CORROSIVE | | * | 95,05,30,001 |
| SODIUM | | CORROSIVE | STRONG REDUCING AGENT | * | |
| Sodium 4(or 5)-methyl-1H-benzotriazolide | 64665-57-2 | CORROSIVE | BASE & REDUCING AGENT | * | 080229001/MSDS |
| SODIUM ACETATE ANHYDROUS | 127-09-3 | IRRITANT | ACID SALT | * | |
| SODIUM ALUMINATE | 1302-42-7 | CORROSIVE | | * | 06.05.04.002/FDS |
| SODIUM AMIDURE | 7782-92-5 | FLAMMABLE / CORROSIVE | | * | 95,11,28,002 |
| SODIUM ASCORBATE | 134-03-2 | | NOT VERY BASIC | * | |
| SODIUM AZIDE | 26628-22-8 | VERY TOXIC | | * | 05,04,18,001/FDS |
| SODIUM BENZOATE | 532-32-1 | HARMFUL | NOT VERY BASIC | * | |
| SODIUM BICARBONATE | 144-55-8 | | SODIUM SALT | * | |
| SODIUM BICHROMATE | 7789-12-0 | TOXIC/IRRITANT/DANGEROUS TO THE ENVIRONMENT | SODIUM SALT | * | |
| SODIUM BIFLUORIDE | 1333-83-1 | TOXIC, CORROSIVE | FLUORIDE | * | 100923001/FDS |
| SODIUM BISULFATE | 10034-88-5 | CORROSIVE | SODIUM SALT | * | |
| SODIUM BISULFITE | 7681-57-4 | CORROSIVE/HARMFUL | SODIUM SALT | * | |
| SODIUM BOROHYDRIDE | 16940-66-2 | SEE SODIUM TETRAHYDROBORATE | | * | |
| SODIUM CARBONATE | 497-19-8 | IRRITANT | | * | 95,05,30,011 |
| SODIUM CHLORATE | 2146053 | HARMFUL/OXIDIZING | SODIUM SALT | * | |
| SODIUM CHLORIDE | 7647-14-5 | | SODIUM SALT | * | 98,08,17,021 |
| SODIUM CHLORITE | 7758-19-2 | VERY TOXIC/CORROSIVE/OXIDIZING/HYGROSCOPIC | SODIUM SALT | * | |
| SODIUM CHLOROACETATE | 3926-62-3 | TOXIC / IRRITANT | | * | 95,05,30,027 |
| SODIUM CHROMATE | 2146108 | TOXIC/DANGEROUS TO THE ENVIRONMENT | OXIDIZER | * | |
| SODIUM DIHYDROGEN PHOSPHATE ANHYDROUS | 7558-80-7 | | | * | |
| SODIUM DIHYDROGEN PHOSPHATE DIHYDRATE | 13472-35-0 | | | * | |
| SODIUM DIHYDROGEN PHOSPHATE MONOHYDRATE | 10049-21-5 | | | * | |
| SODIUM DISULFITE | 7681-57-4 | SEE SODIUM PYROSULFITE | | * | |
| SODIUM DITHIONITE | 7775-14-6 | HARMFUL | REDUCING AGENT | * | |
| SODIUM FLUORIDE | 7681-49-4 | TOXIC | | * | |
| SODIUM HYDRO-2-ETHYL SULFONATE 50 % | | CORROSIVE | CORROSIVE | * | 95,06,01,004 |
| SODIUM HYDROGEN SULFITE | 7631-90-5 | CORROSIVE/HARMFUL | | * | |
| SODIUM HYDROSULFITE | 7775-14-6 | SEE SODIUM DITHIONITE | | * | |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | |
|----------------------------------|------------|--|-----------------------|---|-------------------|
| SODIUM HYDROXIDE | 1310-73-2 | SEE SODA | | | |
| SODIUM HYPOCHLORITE | 7681-52-9 | SEE BLEACH | | | |
| SODIUM HYPOCHLORITE LIQUOR | | CORROSIVE | OXIDIZER | * | 070810003/MSDS |
| SODIUM HYPOPHOSPHITE | 10039-56-2 | | | * | 95,05,30,022 |
| SODIUM HYPOSULFITE | 7775-14-6 | SEE SODIUM DITHIONITE | | | |
| SODIUM ISETHIONATE 50 % | | SEE SODIUM HYDRO-2 ETHYL SULFONATE 50% | | | |
| SODIUM METABISULFITE | 7681-57-4 | SEE SODIUM DISULFITE | | | |
| SODIUM METASILICATE ANHYDROUS | 6834-92-0 | CORROSIVE | | * | |
| SODIUM METASILICATE PENTAHYDRATE | 10213-79-3 | HARMFUL/IRRITANT | | * | |
| SODIUM METHYL DITHIOCARBAMATE | 137-42-8 | HARMFUL | | * | 04,02,06,001/MSDS |
| SODIUM METHYLATE | 124-41-4 | CORROSIVE | | * | 03,02,25,001/FDS |
| SODIUM NITRITE | 7632-00-0 | TOXIC/OXIDIZING | OXIDIZER | * | |
| SODIUM OXALATE | 62-76-0 | CORROSIVE | BASE | * | 01,01,23,002 |
| SODIUM PERBORATE TETRAHYDRATE | 10486-00-7 | HARMFUL | | * | |
| SODIUM PEROXIDE | 1313-60-6 | FLAMMABLE CORROSIVE/OXIDIZING | | * | |
| SODIUM PHOSPHATES | | | SODIUM SALT | * | |
| SODIUM PYROSULFITE | 7681-57-4 | IRRITANT | SODIUM SALT | * | |
| SODIUM SACCHARINATE | | | SODIUM SALT | * | |
| SODIUM SILICATE | 1344-09-8 | CORROSIVE | SODIUM SALT | * | |
| SODIUM SULFATE ANHYDROUS | 7757-82-6 | | SODIUM SALT | * | |
| SODIUM SULFATE DECAHYDRATE | 7727-73-3 | | SODIUM SALT | * | |
| SODIUM SULFITE ANHYDROUS | 7757-83-7 | HARMFUL/IRRITANT | | * | |
| SODIUM SULFITE HYDRATE | 1313-84-4 | CORROSIVE | | * | |
| SODIUM TETRABORATE ANHYDROUS | 1330-43-4 | IRRITANT | SODIUM SALT | * | |
| SODIUM TETRABORATE DECAHYDRATE | 1303-96-4 | | SODIUM SALT | * | |
| SODIUM TETRAHYDROBORATE | 16940-66-2 | TOXIC, CORROSIVE | REDUCING AGENT | * | 100812001/MSDS |
| SODIUM THIOSULFATE ANHYDROUS | 7772-98-7 | IRRITANT/HYGROSCOPIC | | * | |
| SODIUM THIOSULFATE PENTAHYDRATE | 10102-17-7 | IRRITANT | | * | |
| SOUDINOX | | TOXIC, CORROSIVE | MELANGE FLUO-NITRIQUE | * | 071220001-FDS |
| SOYBEAN MILL FEED | | | | * | 04,12,07,010/FDS |
| SPEC TAK | | CORROSIVE | BASE | * | |
| SPECIAL OILS | | | SOLVENTS | * | |
| SPECTRUS BD1512 | | IRRITANT | BIODISPERSANT | * | 061004001/MSDS |
| SPECTRUS NX1164 | | TOXIC/CORROSIVE/PROBABLY SENSITIZER | BIOCIDE | * | 061004002/MSDS |
| SPEZYME AA | | HARMFUL | | * | |
| SPIRAX SARCO BETA 5 | | IRRITANT | CORROSION INHIBITOR | * | 090325001/MSDS |
| SPIRAX SARCO DELTA 986 | | HARMFUL,SENSIBILISANT | BIOCIDE | * | 090327001/MSDS |
| SPIRAX SARCO DELTA 987 | | TOXIC,CORROSIVE | ACID | * | 090325002/MSDS |
| STADIS 425 | | FLAMMABLE,HARMFUL | FUEL ADDITIVE | * | 090204008/MSDS |
| STANNIC CHLORIDE | 7646-78-8 | CORROSIVE | | * | 02,08,30,004 |
| STARCHES | | HARMLESS | HARMLESS | * | |
| STEARIC ACID | 57-11-4 | IRRITANT | ACID | * | 95,05,29,004 |
| STEARYL ALCOHOL | 112-92-5 | SEE 1-OCTADECANOL | | | |
| STYRENE | 100-42-5 | HARMFUL/IRRITANT/FLAMMABLE | | * | |
| SU-8 SERIES RESIST | | IRRITANT | | * | 04,04,06,007/MSDS |
| SUCCINIC ANHYDRIDE | 108-30-5 | IRRITANT | | * | 95,05,30,026 |
| SULFADIAZINE | 68-35-9 | HARMFUL - POSSIBLE SENSITIZER | | * | 04,12,07,007/MSDS |
| SULFAMIC ACID | 5329-14-6 | IRRITANT | WEAK ACID | * | |
| SULFANILIC ACID | 121-57-3 | HARMFUL | | * | 95,06,19,005 |
| SULFOCHROMIC MIX | | CORROSIVE | CORROSIVE | * | 96,05,07,002 |
| SULFUR DIOXIDE | 2025884 | TOXIC/IRRITANT | GAS | * | |
| SULFUR TRIOXIDE | 2025949 | SEE SULFURIC ACID | | | |
| SULFURIC ACID | 8014-97-9 | CORROSIVE | STRONG ACID | * | 90,09,24,001 |
| SULFUROUS ANHYDRIDE | 2025884 | SEE SULFUR DIOXIDE | | | |
| SULFURYL CHLORIDE | 7791-25-5 | CORROSIVE, reacts with water | ACID RELEASE | * | 071213002-MSDS |
| SULPHUR | 7704-34-9 | IRRITANT, SENSITIZER | | * | 1001126001/MSDS |
| SULPHUR HEXAFLUORIDE | 2551-62-4 | TOXIC, IRRITANT | ACID | * | 100422001/MSDS |
| SUPERGREEN (SP98) | 86290-81-5 | FLAMMABLE/TOXIC | HYDROCARBURES MIX | * | 99,01,28,001 |
| SUPRASEC 2025 | | HARMFUL | | * | |
| SUPRASEC 2028 | | HARMFUL | | * | |
| SUPRASEC VM50 REYDEL | | | | * | |
| SURFACE PREPARATING AGENT P50 | | FLAMMABLE, IRRITANT | SOLVENT for PAINT | * | 080326006/FDS |
| SYSTOCHROMAT 1653/1 | | TOXIC/CORROSIVE | | * | 04,03,15,002/FDS |
| TALC | 14807-96-6 | | | * | |
| TALLOW DIMETHYL FAT AMINE | | CORROSIVE / HARMFUL | | * | 95,06,01,014 |
| TANNIC ACID | 72401-53-7 | SEE TANNIN | | | |
| TANNIN | 72401-53-7 | | WEAK ACID | * | |
| TAR | | HARMFUL/FLAMMABLE | NOXIOUS/FLAMMABLE | * | |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | | |
|--|-------------|--|---------------------------|---|---|-------------------|
| TARTRIC ACID | 133-37-9 | IRRITANT | WEAK ACID | * | | |
| TBTU | 125700-67-6 | IRRITANT | | * | | 03,04,01,001/FDS |
| TDI | 584-84-9 | SEE TOLYLENE-2,4-DIISOCYANATE | | | | |
| TDMAT | 3275-24-9 | CORROSIVE/FLAMMABLE | | * | | 02,05,15,002/FDS |
| TEAM "LOW Na + " LIQUID ACID NEUTRALIZER | | IRRITANT | | * | | 04,04,06,009/MSDS |
| TECHNOVIT 2000 LC | | IRRITANT | methacrylates preparation | * | | 02,07,05,002 |
| TECHNOVIT 2000 LC VERNIS | | IRRITANT | methacrylates preparation | * | | 02,07,05,001 |
| TEGO ALKANOL C20P | | | | * | | 03,05,16,013/FDS |
| TEMAZ | 175923-04-3 | see Tetrakis (ethylmethylamido) Zirconium | | * | | 080807001/MSDS |
| TER-BUTYL 2-ETHYL PERHEXANOATE | | | | * | | 92,02,10,002 |
| TERBUTYL PERBENZOATE | 614-45-9 | HARMFUL / FLAMMABLE | | * | | 92,02,10,002 |
| TEREPHTHALIC ALDEHYDE | 623-27-8 | IRRITANT | | * | | 95,05,30,029 |
| TETRAAMINEPLATINUM(II) DINITRATE 1,36% | 20634-12-2 | EXPLOSIVE, IRRITANT, SENSITIZER | | * | | 09/1215001/NEP |
| TETRABUTYLAMMONIUM FLUORIDE (IN THF) | 429-41-4 | CORROSIVE | BASE | * | | 101005001/MSDS |
| TETRACHLOROETHYLENE | 127-18-4 | HARMFUL/DANGEROUS TO THE ENVIRONMENT | | * | | 95,03,23,002 |
| TETRADECANOIC ACID | 544-63-8 | SEE MYRISTIC ACID | | | | |
| TETRAETHYLENEGLYCOL DIACRYLATE | 109-17-1 | HARMFUL / CORROSIVE | | * | | 95,11,28,001 |
| TETRAHYDROFURANNE | 109-99-9 | FLAMMABLE/IRRITANT | SOLVENT | * | | 04,03,15,001/FDS |
| TETRAHYDROTHIOPHENE | 110-01-0 | HARMFUL/IRRITANT/FLAMMABLE | | * | | |
| Tetrakis (ethylmethylamido) Zirconium | 175923-04-3 | REACTS WITH WATER | REDUCING AGENT | * | | 080807001/MSDS |
| TETRAKIS(DIMETHYLAMINO)TITANIUM | | SEE TDMAT | | | | |
| TETRAMETHYL ORTHOSILICATE | 681-84-5 | TOXIC,IRRITANT | | * | | 090827001/NEP |
| TETRAMETHYLAMMONIUM HYDROXIDE | 75-59-2 | TOXIC/CORROSIVE | TOXIC/CORROSIVE | * | | 070327001/NEP |
| TETRAMETHYLENESULFIDE | 110-01-0 | SEE TETRAHYDROTHIOPHENE | | | | |
| TETRAPROPANE | | IRRITANT | SOLVENT | * | | |
| THF | | see TETRAHYDROFURANNE | | | | |
| THIOGLYCOLIC ACID | 68-11-1 | TOXIC / CORROSIVE | | * | | 95,06,19,006 |
| THIOLANE | 110-01-0 | SEE TETRAHYDROTHIOPHENE | | | | |
| THIONYL CHLORIDE | 7719-09-7 | CORROSIVE | | * | | 98,08,17,022 |
| THIOPHANE | 110-01-0 | SEE TETRAHYDROTHIOPHENE | | | | |
| THIOUREE | 62-56-6 | HARMFUL | BASE | * | | |
| THREEBOND 30F-187 | | IRRITANT | | * | | |
| TIN (II) FLUORIDE | 7783-47-3 | HARMFUL, IRRITANT | SOLID | * | * | 081006001/MSDS |
| TINOPAL CBS-X | 27344-41-8 | IRRITANT | | * | | 02,11,28,001/FDS |
| TITAN 2000 EMULSION | | IRRITANT | OXIDISER | * | | 081006001/MSDS |
| TITANIUM TETRACHLORIDE | 7550-45-0 | CORROSIVE | LEWIS ACID | * | | 04,02,17,001/ECH |
| TITANIUM(IV) CHLORIDE | 7550-45-0 | CORROSIVE | LEWIS ACID | * | | 98,04,23,001 |
| TITRIPLEX SOLUTION B | | | | * | | 04,07,01,002 |
| TMAH | 75-59-2 | see TETRAMETHYLAMMONIUM HYDROXIDE | | * | | |
| TMCS | 75-77-4 | SEE TRIMETHYLCHLOROSILANE | | | | |
| TMOS | 681-84-5 | SEE TETRAMETHYL ORTHOSILICATE | | * | | 090827001/NEP |
| TMT 15 | | IRRITANT | | * | | 95,11,29,001 |
| TOLUENE | 108-88-3 | HARMFUL FLAMMABLE | SOLVENT | * | | |
| TOLYLENE-2,4-DIISOCYANATE | 584-84-9 | TOXIC/IRRITANT | | * | | |
| TPH | | SEE TETRAPROPANE | | | | |
| TRANS-2-BUTENOIC ACID | 107-93-7 | SEE CROTONIC ACID | | | | |
| trans-DICHLOROETHYLENE | 156-60-5 | FLAMMABLE, HARMFUL, Dangerous to the environnement | SOLVENT | * | | 040701004-MSDS |
| TRANS-LC TM | 156-60-5 | SEE DICHLOROMETHANE | | | | 04,07,01,004 |
| TRASAR 77395 | | CORROSIVE | BASE | * | | 090819004/FDS |
| TRIALYLAMINE | 102-70-5 | SEE TRI-2-PROPENYLAMINE | | | | |
| TRIBROMOMETHANE | 75-25-2 | SEE BROMOFORME | | | | |
| TRIBUTYL TIN CHLORIDE | 1461-22-9 | | | * | | 05,04,18,002/FDS |
| TRICHLOROACETIC ACID | 76-03-9 | CORROSIVE | STRONG ACID | * | | 98,08,17,015 |
| TRICHLOROETHYLENE | 79-01-6 | HARMFUL | SOLVENT | * | | |
| TRICHLOROMETHANE | 67-66-3 | SEE CHLOROFORME | | | | |
| TRICHLOROMETHYLSILANE | 75-79-6 | FLAMMABLE/REACTS WITH WATER | | * | | 06,08,16,001/EXP |
| TRITHANOLAMINE | 102-71-6 | IRRITANT | AMINE | * | | |
| TRIETHYLALUMINIUM | 97-93-8 | CORROSIVE/FLAMMABLE | LEWIS BASE | * | | 99,07,21,001 |
| TRIETHYLAMINE | 121-44-8 | IRRITANT/HARMFUL/FLAMMABLE | | * | | 95,05,30,023 |
| TRIETHYLAMINE TRISHYDROFLUORID | 73602-61-6 | TOXIC/CORROSIVE | | * | * | 98,12,16,001 |
| TRIETHYLENEGLYCOL | 112-27-6 | IRRITANT | | * | | |
| TRIFLUOROACETIC ACID | 76-05-1 | CORROSIVE | STRONG ACID | * | | |
| TRIGONOX 141 | | | | * | | |
| TRIGONOX 21-C50 | | FLAMMABLE/EXPLOSIVE | | * | | 02,08,29,002 |
| TRIGONOX 36-CD75 | | HARMFUL/IRRITANT | | * | | 02,08,29,001 |
| TRIGONOX B | | | | * | | 02,08,29,003 |
| TRIGONOX C-C75 | | IRRITANT | | * | | 02,08,29,004 |
| TRISOPROPANOLAMINE | 122-20-3 | IRRITANT | | * | | 95,05,29,009 |

TESTED PRODUCT LIST - FIRST EDITION : 15/03/95

| | | | | | | |
|---|------------|-----------------------------------|----------------------------------|---|---|-------------------|
| TRIMERCAPTOTRIAZINE TRISODIUM SALT (15 %) | 17766-26-6 | IRRITANT | | * | | |
| TRIMETHOPRIM | 738-70-5 | TOXIC IF SWALLOWED | | * | | 04,12,07,002/MSDS |
| TRIMETHYLAMINE 40 wt. % solution in water | 75-50-3 | HARMFUL/CORROSIVE/FLAMMABLE | | * | | |
| TRIMETHYLCHINON | 935-92-2 | TOXIC/IRRITANT | | * | | 03,11,04,001/MSDS |
| TRIMETHYLCHLOROSILANE | 75-77-4 | FLAMMABLE CORROSIVE | | * | | 91,01,03,002 |
| TRIPOTASSIUM PHOSPHATE ANHYDROUS | 27176-10-9 | CORROSIVE | | * | | |
| TRIPOTASSIUM PHOSPHATE MONOHYDRATE | 7778-53-2 | CORROSIVE | | * | | |
| TRIPROPYLENEGLYCOL DIACRYLATE | 42978-66-5 | IRRITANT/SENSITIZER | | * | | 06,02,02,001/MSDS |
| TRIS (2-HYDROXYPROPYL)AMINE | 122-20-3 | SEE TRIISOPROPANOLAMINE | | | | |
| TRIS(HYDROXYMETHYL)-AMINOMETHANE | 77-86-1 | IRRITANT | | * | | 95,06,26,013 |
| TRISODIUM PHOSPHATE DODECAHYDATE | 10101-89-0 | CORROSIVE | | * | | |
| TRISODIUM PHOSPHATE HEXAHYDRATE | 15819-50-8 | CORROSIVE | | * | | |
| TRIXYL PHOSPHATE | 25155-23-1 | DANGEROUS TO THE ENVIRONMENT | | * | | 060927001/MSDS |
| TROCLOSENE SODIUM | 2893-78-9 | IRRITANT | STRONG OXIDISER | * | | 071213001-MSDS |
| TRPGDA | | SEE TRIPROPYLENEGLYCOL DIACRYLATE | | | | |
| TURCO 4104 | | CORROSIVE/TOXIC | | * | | 98,08,21,001 |
| TURCO ALUMIGOLD FLUSSIG | | CORROSIVE/TOXIC | | * | | 98,08,21,004 |
| TURCO LIQUID SMUT-GO NC | | CORROSIVE/TOXIC | | * | | 98,08,21,003 |
| TURCO SMUT-GO#4 | | CORROSIVE/TOXIC | | * | | 98,08,21,002 |
| TW1000 RSG-1 OC | | IRRITANT (PEPPER SPRAY) | | * | | 05,06,02,001/MSDS |
| TW1000 RSG-2 OC | | IRRITANT (PEPPER SPRAY) | | * | | 05,06,02,001/MSDS |
| TW1000 RSG-4 OC | | IRRITANT (PEPPER SPRAY) | | * | | 05,06,02,001/MSDS |
| TW1000 RSG-5 OC | | IRRITANT (PEPPER SPRAY) | | * | | 05,06,02,001/MSDS |
| TYPHON PREMIUM | | IRRITANT | PAINT | * | | 080326001/FDS |
| UFACID KW | | CORROSIVE | | * | | 02,02,19,005 |
| UGILEC T | | HARMFUL | | * | | |
| UNICHROME CR 842 | | VERY TOXIC/CMR/CORROSIVE | | * | * | 061023002/MSDS |
| UNLEADED FUEL | 86290-81-5 | SEE SUPERGREEN (SP98) | | * | | 99,02,22,001 |
| UREA | 57-13-6 | | BASE | * | | |
| VEGETABLE COAL | | | | * | | 95,06,26,006 |
| VERNIS 10-140 | | | RESIN | * | | 89,09,04,003 |
| VIAPAL UP 002/60 | | | POLYESTER | * | | 05,07,07,001/MSDS |
| VIKLOX R 18 EPOXIDE | | TOXIC/IRRITANT | | * | | 95,06,01,008 |
| VINYL ACETATE MONOMER | 108-05-4 | HARMFUL | | * | | 95,05,30,019 |
| VINYL PARATERTIOBUTYL BENZOATE | | HARMFUL / IRRITANT | | * | | 95,05,30,032 |
| VITEC 3000 | | IRRITANT | | * | | |
| VITEC 8000 | 22042-96-2 | IRRITANT | | * | | 04,07,01,002 |
| VITRIDE | | | STRONG REDUCING AGENT | * | | 95,06,02,002 |
| VOLPO C20 | 9004-95-9 | HARMFUL | | * | | 03,05,16,012/FDS |
| WBC 34 SPEZIALLÖSER | | IRRITANT | SOLVENT | * | | 01,01,23,005 |
| WELD-GUARD STAINLESS STEEL PICKLING GEL | | TOXIC, CORROSIVE | MELANGE FLUO-NITRIQUE | * | * | 071217001-MSDS |
| WET SODIUM SALT CAKE | | CORROSIVE/TOXIC | | * | | 01,01,23,001 |
| WHITE-SPIRIT | | FLAMMABLE | 17% AROMATICS | * | | 95,03,23,001 |
| XP SU-8 DEVELOPER | | | | * | | 04,06,24,001/MSDS |
| XP SU-8 SERIES RESISTS | | | | * | | 04,06,24,001/MSDS |
| XYLENES | | HARMFUL/IRRITANT/FLAMMABLE | | * | | |
| ZESTRON SD 300 | | FLAMMABLE | | * | | 04,03,12,001/MSDS |
| ZETAG 8818 | | IRRITANT | POLYMERES | * | | 090715001/FDS |
| ZETAG 8868FS | | IRRITANT | POLYMERES | * | | 090715002/FDS |
| ZINC CHLORIDE | 7646-85-7 | CORROSIVE | ZINC SALT | * | | |
| ZINC DIETHYLE | 557-20-0 | CORROSIVE, REACTS WITH WATER | REDUCING AGENT, ALKYLATING AGENT | * | | 060602001/MSDS |
| ZINC SULFATE HEPTAHYDRATE | 7446-20-0 | IRRITANT | | * | | 95,06,26,012 |
| ZINC SULFATE MONOHYDRATE | 7446-19-7 | IRRITANT | | * | | 95,06,26,012 |
| ZINEB | 12122-67-7 | IRRITANT, SENSITIZER | PESTICIDE | * | | 100412008/MSDS |